

HOSPITAL PRODUCT CATALOG **2026**

NYOrtho

Transforming
PATIENT CARE
Through Innovative Design and Manufacturing

Find out why some of the world's largest healthcare companies private label by NYOrtho....

QUALITY



Commitment to
Excellence



Quality
Control



Product
Development



Manufacturing
Expertise

FLEXIBILITY



Customer Specific
Product



Collaboration



Customer Care



Short Runs

DEPENDABILITY



Responsiveness



Shorter Lead Time



Capability



Accuracy

Our **client-first attitude**
sets us apart because we possess the capability to turn your
concept into reality within a good time frame.

PROVIDING PATIENT SAFETY & COMFORT

THROUGH
DESIGN &
INNOVATION



PRODUCT CODING

Look for the informative icons for product enhancements.

VF

VISCO MEMORY FOAM

AM

ANTI-MICROBIAL

WR

WATER-RESISTANT

CSA

CUSTOM SIZING AVAILABLE

HA

HYPOALLERGENIC

WP

WATER-PROOF

TC

THERMAL CONDUCTIVITY

WOUND MANAGEMENT

pages 4-13

SEATING & POSITIONING

pages 14 - 29

WHEELCHAIR & GERI-CHAIR ACCESSORIES

pages 30 - 37

BED SAFETY & ACCESSORIES

pages 38 - 46

THERAPY & AMBULATION

pages 47 - 50

PATIENT SAFETY

pages 51

ORTHOPEDIC REHABILITATION

pages 52 - 55

OCCUPATIONAL & ERGONOMIC

pages 56 - 57

PROTECTIVE SAFETY WEAR

pages 58 - 60

CONTACT US

800.699.2358
cs@nyorthousa.com

ANTI-MICROBIAL

ZERO-G SERIES

RECOMMENDED USE: Individuals unable to avoid long periods of uninterrupted pressure over bony prominences are at increased risk for the development of necrosis and ulceration.

ZERO-G Boot™

ZERO-G Heel Cloud™

ZERO-G PERFORMA BOOT™

ZERO-G PERFORMA GTO BOOT™

ZERO-G Heel Pillow™

ZERO-G Heel Pillow ULTRA™

ZERO-G Heel Cloud Cushion™



ZERO-G Boot™ AM

The **NYOrtho ZERO-G Boot™** combines a cutting-edge design with high-performance fabrics to eliminate pressure and friction to the skin. There are many lookalikes but only one brand features an **ANTIMICROBIAL FABRIC**. We've made our product stand out from the pack and changed the color of our popular Zero-G Boot so you won't make a mistake and buy an inferior product.

MODEL	SIZE	FITS	WIDTH	UOM
9518-S	Small	Petite Adult / Pediatric	5-10 in.	Each
9518-M	Medium	Average Adult	11-20 in.	Each
9518-L	Large	Large Adult / Bariatric	21-25 in.	Each
FA-GEL-4x8	Replacement Gel-Pack	To choose the correct size, measure the circumference of the lower part of the calf muscle.		
9518-BAG	Washing Machine Bag			



Independent laboratory tests via C-PicoPress Interface Pressure Monitor, with ambient temperature consistent with ISO 554-1976E, prove pressure is effectively removed from the heel.

EC Service Project#1385

SUPER-SOFT VELVET LINING

SUPER-SOFT VELVET cloth protects thin fragile skin from tears caused by trauma, reducing friction and shear.

HIGH PERFORMANCE FABRIC

Antimicrobial, mildew resistant & water-resistant.

ADAPTABLE STRAP

Permits this device to be secured at two alternate positions. Lower position controls dorsiflexion.

PADDED HOOK AND LOOP

Easy-on, easy-off padded closure allows easy application and removal with a comfortable fit.

FOAM CUSHIONED SIDES

Lateral and medial malleoli – bony prominence on both sides of the ankle, are softly cushioned.

ZERO PRESSURE

Extended, cushioned design floats the heel over an air cavity for **ZERO** pressure and provides air circulation to aid the healing of existing ulcers.

A REMOVABLE GEL-PACK

May be utilized to further lower pressure beneath the ankle and keeps feet cool. (FA-GEL-4x8)



ZERO-G Performa Boot™



The NYOrtho ZERO-G Performa Boot™ combines an aggressive design with high-performance fabrics to eliminate pressure and friction to the skin. **NEW PERFORMANCE COVER + ZERO-G Boot!**

MODEL	SIZE	FITS	WIDTH	UOM
9518P-S	Small	Petite Adult / Pediatric	5-10 in.	Each
9518P-M	Medium	Average Adult	11-20 in.	Each
9518P-L	Large	Large Adult / Bariatric	21-25 in.	Each
FA-GEL-4x8	Replacement Gel-Pack	To choose the correct size, measure the circumference of the lower part of the calf muscle.		
9518-BAG	Washing Machine Bag			



WIPE CLEAN PERFORMANCE COVER: W/ADDED ANTIMICROBIAL, MILDEW RESISTANT, WATER RESISTANT PROTECTION

High performance, pressure reducing wipe clean fabric, promotes healing via active evaporation of moisture and cooling at patient pressure points.



ZERO-G Performa GTO Boot™



The NYOrtho ZERO-G Performa GTO Boot™, improves on the wipe clean ZERO-G Performa Boot with a new **EASY-TO-ACCESS-PANEL** that simplifies heel alignment, promotes skin assessment, facilitates application of medicine and eliminates discomfort.

MODEL	SIZE	FITS	WIDTH	UOM
9518P-GTO-S	Small	Petite Adult / Pediatric	5-10 in.	Each
9518-GTO-M	Medium	Average Adult	11-20 in.	Each
9518-GTO-L	Large	Large Adult / Bariatric	21-25 in.	Each
FA-GEL-4x8	Replacement Gel-Pack	To choose the correct size, measure the circumference of the lower part of the calf muscle.		
9518-BAG	Washing Machine Bag			



ELIMINATES DISCOMFORT

Easy to access door eliminates patient discomfort when removing and reapplying the boot that may irritate sensitive skin.

SKIN ASSESSMENT

Simplifies skin assessment and the application of medicine without the need to remove the boot.

HEEL ALIGNMENT

EASY-TO-ACCESS panel makes it easy to align the heel over the air cavity when fitting the patient.



ZERO-G Heel Cloud™

AM

WR

US Patent 11,833,070

The ZERO-G Heel Cloud™ combines craftsmanship and innovation providing a sides-down, open door **patented** design that provides an accessible, easy to monitor wound site. High-performance fabrics promote and maintain a cool and clean healing surface. Hypoallergenic micro-fiber filling, comfortably cradles the foot, eliminating pressure and friction to the skin while offloading the heel and bony prominence on both sides of the ankle.

MODEL	DESCRIPTION	SIZE	UOM
9538	Standard	Universal	Each
9538W	W/Optional Wedge	Universal	Each

THREE DVT HOLES

Left, right and center for attachment of DVT pressure devices.



WIPE CLEAN NYLON

Exterior fabric is easy to maintain.

BREATHABLE, WATER-RESISTANT, ANTIMICROBIAL

Peach-skin fabric helps promote air circulation and wick moisture away.

SOFT INTERIOR PEACHSKIN

Fabric with texture gently helps maintain heel position minimizing patient movement.

HYPOALLERGENIC MICRO-FIBER FILLING

Comfortably cradles the foot and cushions the ankles, eliminating pressure while offloading the heel.

OPTIONAL WEDGE

Stabilizer helps prevent foot and leg rotation.
MODEL 9538W





OPEN BOTTOM OFFLOADS THE HEEL
Eliminating pressure and thereby pressure sores.



SIDES-DOWN, OPEN DOOR DESIGN
Provides an accessible, easy to monitor wound site and prevents plantar flexion with an easy to adjust positioner.

ELIMINATE PLANTAR FLEXION

Prevents plantar flexion with an easy to adjust positioner.






ZERO-G Heel Pillow Ultra™



The **NYOrtho ZERO-G Heel Pillow ULTRA™** transforms everyday therapy providing exceptional support and comfort with a new contoured design that conforms to the leg's natural shape and position, while off-loading the heel into the air. The **ULTRA** provides a more relaxing and restful experience promoting patient position, while preventing, and providing relief from pressure sores.

Meets CA 117 standards

MODEL	STYLE	SIZE	UOM
9901	Single Leg	12" x 17" x 4¾" x 3½"	Each
9902	Both Legs	24" x 17" x 4¾" x 3½"	Each



PRESERVES PATIENT POSITION

New ULTRA design preserves patient position, while relaxing the spine and decreasing lower back pain.

CONVOLUTED FOAM

Convoluted foam topper increases ventilation, optimizing pressure relief and comfort.

OFF LOADS HEELS

Cushion design off-loads heels into the air, preventing and providing relief from pressure sores, while securely and comfortably positioning both legs.

HIGH-PERFORMANCE, WIPE CLEAN COVER

Breathable, low-shear, stretch cover with **ANTIMICROBIAL** protection reduces friction. Slip resistant bottom prevents sliding.

BI-DIRECTIONAL

Both styles can be positioned in either direction to maintain slightly lifted or lowered position of the heel while providing abduction.





ZERO-G Heel Pillow™ AM WR

The NYOrtho ZERO-G Heel Pillow™ is specially designed to off-load the heels into the air, relaxing the spine and decreasing lower back pain, thereby preventing, and providing relief from pressure sores.

MODEL	STYLE	SIZE	UOM
9900	Standard	24" x 16" x 7" x 5½"	Each
9900S	With straps	24" x 16" x 7" x 5½"	Each
9900-G	Gel	24" x 16" x 7" x 5½"	Each

Meets CA 117 standards

OFF LOADS HEELS

Cushion design off-loads heels into the air, preventing and providing relief from pressure sores, while securely and comfortably positioning both legs.

HIGH-PERFORMANCE, WIPE CLEAN COVER

Breathable, low-shear, stretch cover with **ANTIMICROBIAL** protection reduces friction. Slip resistant bottom prevents sliding.

POPLITEAL SUPPORT

Popliteal support design fills in the gap behind the knees for superior comfort, elevation and support.



Heel Cloud™ Cushion AM WR VF

Unlike standard therapeutic cushions, that only provide the lift required to suspend the heel in the air, the **NYOrtho Heel Cloud™ Cushion** exceeds current conventional therapeutic cushions and pillows by providing exceptional cloud-like support and comfort with a new design that conforms to the body's natural position while off-loading the heel into the air. A wide landing pad allows for varied positions thereby increasing patient comfort.

MODEL	SIZE	UOM
9903	30" x 23" x 6"	Each

Meets CA 117-2013 standards

OFF LOADS HEELS

Cushion design off-loads heels into the air, preventing and providing relief from pressure sores, while securely and comfortably positioning both legs.

HIGH-PERFORMANCE, WIPE CLEAN COVER

Breathable, low-shear, stretch cover with **ANTIMICROBIAL** protection reduces friction. Slip resistant bottom prevents sliding.

CONFORMS TO YOUR NATURAL SHAPE

Three layers of high resiliency foam and **VISCOELASTIC** foam conforms to your natural shape.

POPLITEAL SUPPORT

Popliteal support design fills in the gap behind the knees for superior comfort, elevation and support.



HEEL & ELBOW PROTECTORS

RECOMMENDED USE: Patients at risk for skin breakdown and friction burns.

NYOrtho Quilted Heel and Elbow Protectors provide comfort and protection to sensitive skin, reducing pressure and increasing blood flow to the affected area. A soft and fluffy spun fiber-fill helps provide extra comfort and protection. Padded hook & loop closure allows for easy application and removal.

Super Quilted Heel Protector

AM WR

MODEL	W/GEL-PACK	SIZE	DIMENSIONS	UOM
9514	9514-GP	Standard	10" x 10"	Pair
9514L	N/A	Bariatric	12" x 12"	Pair

SUPER FLUFFY FILL

Designed with SUPER FLUFFY, soft, hypoallergenic, micro-fiber filling for added comfort.

HIGH PERFORMANCE FABRIC

Breathable, water-resistant, antimicrobial peach-skin fabric helps promote air circulation and wick moisture away.

PADDED HOOK AND LOOP CLOSURE

For superlative comfort and easy application or removal.

MACHINE WASHABLE



Quilted Heel & Calf Protector

AM WR



MODEL	SIZE	DIMENSIONS	UOM
9515	Universal	14" x 9"	Pair

HIGH PERFORMANCE FABRIC

Breathable, water-resistant, antimicrobial peach-skin fabric helps promote air circulation and wick moisture away.

PADDED HOOK AND LOOP CLOSURE

For superlative comfort and easy application or removal.

MACHINE WASHABLE



Quilted Elbow Protector



MODEL	SIZE	UOM
9360	Universal	Pair

HIGH PERFORMANCE FABRIC

Breathable, water-resistant, antimicrobial peach-skin fabric helps promote air circulation and wick moisture away.

MACHINE WASHABLE.



Ultra Quilted Heel Protector



MODEL	SIZE	DIMENSIONS	UOM
9516	Standard	9" x 10"	Pair
9516L	Bariatric	11" x 11½"	Pair

ULTRA-SOFT

Anti-microbial velvet fabric is a gentle on the skin.

ANTIMICROBIAL

No odor no bacteria: Thanks to the built-in hygiene function



Fun Quilted Heel Protector

MODEL	SIZE	DIMENSIONS	UOM
9511	Standard	9" x 10"	Pair
9511L	Bariatric	11" x 11½"	Pair

COMFORTABLE

Designed with fun cotton fabric – a soft, breathable and durable fabric with excellent wickability.

Color and/or design of fabric subject to change.

ANTI-MICROBIAL

SKIN PROTECTION

RECOMMENDED USE: Patients with fragile skin; to protect IV / wound sites and to protect skin from tears and abrasions.

NYOrtho
Derma Shield™

NYOrtho Derma Shield™ **AM**

Derma Shield™ protective arm and leg sleeves are a non-compression skin covering that help provide a protective barrier against tears and abrasions for fragile skin tissue and adapts to different body shapes. Protective arm and leg sleeves can also help prevent patient disruption of medical treatments by camouflaging both IV sites and wound sites by blending in with the patient's skin tone.



ANTIMICROBIAL



**NOT MADE
FROM NATURAL
RUBBER LATEX**



**MACHINE
WASHABLE**



**PROTECTS
SKIN**



**4-WAY STRETCH
SOFT MATERIAL**



ARM SIZING BANDS



Small



Medium



Large



X-Large

MODEL	STYLE	SIZE	LENGTH	CIRC.	UOM
9402-S	Arm	Small	15½"	8½"	Pair
9402-M	Arm	Medium	16½"	11"	Pair
9402-L	Arm	Large	18"	11"	Pair
9402-XL	Arm	XL	20"	13"	Pair
9100-U	Leg	Universal	13½"	11"	Pair



Economy Sheepskin Boot

RECOMMENDED USE: Patients at low risk of skin breakdown and friction burns.

Synthetic sheepskin products help make everything cozy and comfortable. Soft and plush material feels wonderful against the skin. High-pile, low shear, non-allergenic synthetic sheepskin fur, promotes air circulation and reduces risk of pressure sores.



Other Soft Sheepskin Products



9517

- Boot covers heel, ankle and leg up to mid-calf.
- Non-allergenic, flame retardant, imitation sheepskin.
- Comfortably secured with soft loop straps.
- One size fits all.
- Meets CA 117 standards.



Page 37



Page 39



Page 42

Economy Heel Protector



9509-A



9509-B

Economy Elbow Protector

MODEL	DESCRIPTION	SIZE	UOM
9509-A	Heel Protector	Universal	Pair
9509-B	Elbow Protector	Universal	Pair
9517	Sheepskin Boot	Universal	Pair

ANTI-MICROBIAL & WATER-RESISTANT

WHEELCHAIR CUSHIONS

NYOrtho™ Series
APEX™ Series
APEX CORE™ Series

RECOMMENDED USE: Patients at risk for pressure ulcers.

Three Weight Distribution Designs / Two High Performance Covers

Choose from three different foam/gel distribution plans! The **NYOrtho 9595** series provides a tried and true wheelchair gel-foam cushion design that provides support and weight distribution. The innovative **APEX 9597** series cools and reduces body interface pressure. Combine a firm core foundation with the superlative APEX design and you have the **APEX CORE 9599** series.

NYOrtho™ Series

NYOrtho™ Design Cools and Reduces Body Interface Pressure

Dual density foam layers provide a medium density top entry layer and a higher density base layer for support and pressure distribution.

Aqueous Chambered Gel-Pod

Reduces gel migration and dissipates metabolic heat, continuously maintaining a cool interface between the body and cushion.

APEX™ Series

Innovative Design Cools and Reduces Body Interface Pressure

Dual Density Foam layers provides a soft top entry layer promoting direct interface with cooling gel-pod and a higher density base layer for support and pressure distribution. This top-design brings the gel-pod up close so you can feel the superior comfort and stability.

Aqueous Chambered Gel-Pod

Available dual and quad chambered gel-pod reduces gel migration and dissipates metabolic heat, continuously maintaining a cool interface between the body and cushion.

APEX CORE™ Series

Firm Core Foundation

High Density beveled polyethylene provides a firm foundation that eliminates hammocking.

Innovative Design Cools And Reduces Body Interface Pressure

Dual Density Foam layers provides a soft top entry layer promoting direct interface with cooling gel-pod and a higher density base layer for support and pressure distribution. This top-design brings the gel-pod up close so you can feel the superior comfort and stability.

Aqueous Chambered Gel-Pod

Dual chambered gel-pod reduces gel migration and dissipates metabolic heat, continuously maintaining a cool interface between the body and cushion.

Hi-Performance Wheelchair Cushion Covers



RECOMMENDED USE: Patients at risk for pressure ulcers.
Meets CA 117 standards

MODELS 9595, 9597, 9599

BREATHABLE STRETCH REMOVABLE COVER

High performance pressure reducing fabric promotes healing via active evaporation of moisture and cooling at patient pressure points.

NON-SLIP BOTTOM


Aggressive non-slip bottom with easy to attach quick-release straps keep cushion in place.



**ANTI-MICROBIAL, MILDEW RESISTANT,
FLAME RETARDANT & WATER-RESISTANT**

EASY-TO-ATTACH QUICK-RELEASE STRAPS

Conveniently attaches cushion to wheelchair frame.



Sure-Chek™ Fabric Wheelchair Cushion Covers



RECOMMENDED USE: Patients at risk for pressure ulcers.

MODEL 93-95-XXX-XXXX NYOrtho Series

MODEL 93-97-XXX-XXXX APEX Series

MODEL 93-99-XXX-XXXX APEX CORE Series



• FLAME RESISTANT

All Sure-Chek healthcare fabrics are flame resistant and self-extinguish when a source flame is removed.

• THERMAL CONDUCTIVITY

Sure-Chek healthcare fabrics all perform within a desirable range for patient comfort without sacrificing important technical performance benefits.

• HYPOALLERGENIC

Sure-Chek healthcare fabrics have been engineered to be hypoallergenic and comfortable for the user.

• ANTIMICROBIAL FABRIC PROTECTION

The antimicrobial agent helps protect the fabric from microbial attack and reduces undesirable organic odors.

• FLUID PROOF & MOISTURE VAPOR TRANSMISSION

All Sure-Chek healthcare fabrics are fluid and moisture proof, reducing the risk of contamination of the mattress from liquids or bodily fluids.



Do we manufacture custom sizes? WE CERTAINLY DO!

We can manufacture any kind and size cushion you need. We understand your customers' needs. Doctors, nurses and therapists are always tinkering with foam to make the necessary adjustments to a cushion to fit the specific needs of each patient. They often wish that the manufacturer could make them a cushion just the way they want it.

Our designers are ready to help. There are many different directions a custom product can take. Different covers can be applied to different cushions. Inner foam densities can be selected to your specifications and shapes can be modified to provide an improved product that provides the optimum rehabilitation objective for your patient. Anything you need, we will do it with minimal lead time.



NYOrtho™ Series

Meets CA 117 standards



NYOrtho Gel-Foam Cushions combine a dual chambered gel insert within a dual density foam body. This combination helps prevent pressure injuries, improve patient comfort and correct patient positioning. Features a high-performance stretch cover that improves moisture evaporation with an aggressive non-slip bottom keeping the cushion firmly in place. UOM - each, unless otherwise stated.

RECOMMENDED USE: Patients at risk for pressure ulcers.



NYOrtho Gel-Foam Cushion

Size	2" Height	3" Height	4" Height
16x16 in.	9595-GEL-161602	9595-GEL-161603	9595-GEL-161604
18x16 in.	9595-GEL-181602	9595-GEL-181603	9595-GEL-181604
18x18 in.	9595-GEL-181802	9595-GEL-181803	9595-GEL-181804
20x16 in.	9595-GEL-201602	9595-GEL-201603	9595-GEL-201604
20x18 in.	9595-GEL-201802	9595-GEL-201803	9595-GEL-201804
22x18 in.	9595-GEL-221802	9595-GEL-221803	9595-GEL-221804
24x18 in.	9595-GEL-241802	9595-GEL-241803	9595-GEL-241804
Custom sizes available - any length, width or height.			

NYOrtho Foam Cushion

2" Height	3" Height	4" Height
9595-161602	9595-161603	9595-161604
9595-181602	9595-181603	9595-181604
9595-181802	9595-181803	9595-181804
9595-201602	9595-201603	9595-201604
9595-201802	9595-201803	9595-201804
9595-221802	9595-221803	9595-221804
9595-241802	9595-241803	9595-241804
Custom sizes available - any length, width or height.		

NYOrtho BARIATRIC Gel-Foam Cushion

Size	3" Height	4" Height
26 x 18 in.	9595B-GEL-261803	9595B-GEL-261804
26 x 20 in.	9595B-GEL-262003	9595B-GEL-262004
28 x 18 in.	9595B-GEL-281803	9595B-GEL-281804
28 x 20 in.	9595B-GEL-282003	9595B-GEL-282004
30 x 18 in.	9595B-GEL-301803	9595B-GEL-301804
30 x 20 in.	9595B-GEL-302003	9595B-GEL-302004
32 x 18 in.	9595B-GEL-321803	9595B-GEL-321804
32 x 20 in.	9595B-GEL-322003	9595B-GEL-322004
Custom sizes available - any length, width or height.		

NYOrtho BARIATRIC Foam Cushion


Size	3 Height	4 Height
26 x 18 in.	9595B-261803	9595B-261804
26 x 20 in.	9595B-262003	9595B-262004
28 x 18 in.	9595B-281803	9595B-281804
28 x 20 in.	9595B-282003	9595B-282004
30 x 18 in.	9595B-301803	9595B-301804
30 x 20 in.	9595B-302003	9595B-302004
32 x 18 in.	9595B-321803	9595B-321804
32 x 20 in.	9595B-322003	9595B-322004
Custom sizes available - any length, width or height.		

NYOrtho Convoluted Gel-Foam Cushion



Size	2" Height	3" Height
16x16 in.	9595-CGL-161602	9595-CGL-161603
18x16 in.	9595-CGL-181602	9595-CGL-181603
18x18 in.	9595-CGL-181802	9595-CGL-181803
20x16 in.	9595-CGL-201602	9595-CGL-201603
20x18 in.	9595-CGL-201802	9595-CGL-201803
22x18 in.	9595-CGL-221802	9595-CGL-221803
24x18 in.	9595-CGL-241802	9595-CGL-241803
Custom sizes available - any length, width or height.		

NYOrtho Convoluted Foam Cushion



Size	Foam
16x16x3 in.	9594-161603
18x16x3 in.	9594-181603
20x16x3 in.	9594-201603
Custom sizes available - any length, width or height.	

NYOrtho Wedge Gel-Foam Cushion Wedge Foam Cushion


RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring posture support, that are inclined to slide forward.



Size	Gel-Foam	Foam
16x16x4x2 in.	9595-WGL-161642	9595-W-161642
18x16x4x2 in.	9595-WGL-181642	9595-W-181642
18x18x4x2 in.	9595-WGL-181842	9595-W-181842
20x16x4x2 in.	9595-WGL-201642	9595-W-201642
20x18x4x2 in.	9595-WGL-201842	9595-W-201842
22x18x4x2 in.	9595-WGL-221842	9595-W-221842
24x18x4x2 in.	9595-WGL-241842	9595-W-241842
Custom sizes available - any length, width or height.		

NYOrtho Pommel Gel-Foam Cushion Pommel Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring hip positioning.



Gel-Foam

Size	2" Height	3" Height
16x16 in.	9595-GP-161602	9595-GP-161603
18x16 in.	9595-GP-181602	9595-GP-181603
18x18 in.	9595-GP-181802	9595-GP-181803
20x16 in.	9595-GP-201602	9595-GP-201603
20x18 in.	9595-GP-201802	9595-GP-201803
22x18 in.	9595-GP-221802	9595-GP-221803
24x18 in.	9595-GP-241802	9595-GP-241803
Custom sizes available - any length, width or height.		

Foam

Size	2" Height	3" Height
16x16 in.	9595-P-161602	9595-P-161603
18x16 in.	9595-P-181602	9595-P-181603
18x18 in.	9595-P-181802	9595-P-181803
20x16 in.	9595-P-201602	9595-P-201603
20x18 in.	9595-P-201802	9595-P-201803
22x18 in.	9595-P-221802	9595-P-221803
24x18 in.	9595-P-241802	9595-P-241803
Custom sizes available - any length, width or height.		

NYOrtho Wedge Pommel Gel-Foam Cushion Wedge Pommel Foam Cushion


RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring hip positioning.



Size	Gel-Foam	Foam
16x16x4x2 in.	9595-GWP-161642	9595-WP-161642
18x16x4x2 in.	9595-GWP-181642	9595-WP-181642
18x18x4x2 in.	9595-GWP-181842	9595-WP-181842
20x16x4x2 in.	9595-GWP-201642	9595-WP-201642
20x18x4x2 in.	9595-GWP-201842	9595-WP-201842
22x18x4x2 in.	9595-GWP-221842	9595-WP-221842
24x18x4x2 in.	9595-GWP-241842	9595-WP-241842
Custom sizes available - any length, width or height.		

NYOrtho Coccyx Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers in the coccyx are.



Size	3" Height
16x16x3 in.	9595-C-161603
18x16x3 in.	9595-C-181603
18x18x3 in.	9595-C-181803
20x16x3 in.	9595-C-201603
20x18x3 in.	9595-C-201803
22x18x3 in.	9595-C-221803
24x18x3 in.	9595-C-241803
Custom sizes available - any length, width or height.	

Custom Cushions Upon Request

Any length, width or height.

NYOrtho Independent Cell Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers. Independent cell foam adapts and provides pressure relief where needed.



Size	3" Height
16x16x3 in.	9595-IC-161603
18x16x3 in.	9595-IC-181603
18x18x3 in.	9595-IC-181803
20x16x3 in.	9595-IC-201603
20x18x3 in.	9595-IC-201803
22x18x3 in.	9595-IC-221803
24x18x3 in.	9595-IC-241803

Solid Seat Insert

RECOMMENDED USE: Patients needing a firm base that eliminates hammocking associated with wheelchairs.



A must have accessory for any wheel chair cushion, The **NYOrtho Solid Seat Insert** provides a firm base that eliminates hammocking associated with wheelchairs. Turns any cushion into a non-slip cushion.

- High density beveled polyethylene provides a firm foundation that eliminates hammocking.
- Aggressive non-slip cover keeps foundation in place.
- Wipe-clean surface.

Size	Style Wedge	Model	UOM
16x16 in.	Full Seat	9-SSI-1616	Each
18x16 in.	Full Seat	9-SSI-1816	Each
18x18 in.	Full Seat	9-SSI-1818	Each
20x16 in.	Full Seat	9-SSI-2016	Each
20x18 in.	Full Seat	9-SSI-2018	Each
22x18 in.	Full Seat	9-SSI-2218	Each
24x18 in.	Full Seat	9-SSI-2418	Each
Custom sizes available - any length, width or height.			

Half Seat Wedge

RECOMMENDED USE: For patients suffering from lateral pelvic tilt requiring a wedge to equalize the pelvis.



The **NYOrtho Half Seat Wedge** provides the support necessary to help correct an uneven pelvis. With the thicker end of the wedge under the patient's lower pelvis, the sitting position is raised to a neutral position,

- Corrects posture of those patients suffering from lateral pelvic tilt.
- High density beveled polyethylene provides a firm foundation to offset lateral pelvic tilt, resulting in increased trunk elongation, upright sitting position and decreased leaning to one side.
- Aids feeding and other activities where posture and alertness are important.
- Aggressive non-slip cover keeps wedge in place.
- Sold as a set of two.

Size	Style Wedge	Model	UOM
9x16 in.	Half Seat	9-HSI-0916	Set
Custom sizes available - any length, width or height.			



APEX™ Series

Meets CA 117 standards



The NYOrtho APEX™ Gel-Foam Cushion features an innovative design that chambers the dual gel-pod closer to the top of the cushion providing a cooler interface between the patient and wheelchair. A high-performance, stretch cover improves moisture evaporation, enhancing APEX™ design features.

RECOMMENDED USE: Patients at risk for pressure ulcers.



APEX™ Gel-Foam Cushion

Size	2" Height	3" Height	4" Height
16x16 in.	9597-GEL-161602	9597-GEL-161603	9597-GEL-161604
18x16 in.	9597-GEL-181602	9597-GEL-181603	9597-GEL-181604
18x18 in.	9597-GEL-181802	9597-GEL-181803	9597-GEL-181804
20x16 in.	9597-GEL-201602	9597-GEL-201603	9597-GEL-201604
20x18 in.	9597-GEL-201802	9597-GEL-201803	9597-GEL-201804
22x18 in.	9597-GEL-221802	9597-GEL-221803	9597-GEL-221804
24x18 in.	9597-GEL-241802	9597-GEL-241803	9597-GEL-241804


Custom sizes available - any length, width or height.

APEX™ Foam Cushion

2" Height	3" Height	4" Height
9597-161602	9597-161603	9597-161604
9597-181602	9597-181603	9597-181604
9597-181802	9597-181803	9597-181804
9597-201602	9597-201603	9597-201604
9597-201802	9597-201803	9597-201804
9597-221802	9597-221803	9597-221804
9597-241802	9597-241803	9597-241804

Custom sizes available - any length, width or height.

APEX QUAD Gel-Foam Cushion



Size	2 in. Height	3 in. Height	4 in. Height
16x16 in.	9597-4GEL-161602	9597-4GEL-161603	9597-4GEL-161604
18x16 in.	9597-4GEL-181602	9597-4GEL-181603	9597-4GEL-181604
18x18 in.	9597-4GEL-181802	9597-4GEL-181803	9597-4GEL-181804
20x16 in.	9597-4GEL-201602	9597-4GEL-201603	9597-4GEL-201604
20x18 in.	9597-4GEL-201802	9597-4GEL-201803	9597-4GEL-201804
22x18 in.	9597-4GEL-221802	9597-4GEL-221803	9597-4GEL-221804
24x18 in.	9597-4GEL-241802	9597-4GEL-241803	9597-4GEL-241804

Custom sizes available - any length, width or height.

APEX Wedge Gel-Foam Cushion Wedge Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring posture support, that are inclined to slide forward.



Size	Gel-Foam	Foam
16x16x4x2 in.	9597-WGL-161642	9597-W-161642
18x16x4x2 in.	9597-WGL-181642	9597-W-181642
18x18x4x2 in.	9597-WGL-181842	9597-W-181842
20x16x4x2 in.	9597-WGL-201642	9597-W-201642
20x18x4x2 in.	9597-WGL-201842	9597-W-201842
22x18x4x2 in.	9597-WGL-221842	9597-W-221842
24x18x4x2 in.	9597-WGL-241842	9597-W-241842

Custom sizes available - any length, width or height.

APEX Pommel Gel-Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring hip positioning.



Size	2" Height	3" Height
16x16x3 in.	9597-GP-161602	9597-GP-161603
18x16x3 in.	9597-GP-181602	9597-GP-181603
18x18x3 in.	9597-GP-181802	9597-GP-181803
20x16x3 in.	9597-GP-201602	9597-GP-201603
20x18x3 in.	9597-GP-201802	9597-GP-201803
22x18x3 in.	9597-GP-221802	9597-GP-221803
24x18x3 in.	9597-GP-241802	9597-GP-241803
Custom sizes available - any length, width or height.		

APEX Coccyx Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers in the coccyx area.



Size	2" Height	3" Height	4" Height
16x16 in.	9597-GC-161602	9597-GC-161603	9597-GC-161604
18x16 in.	9597-GC-181602	9597-GC-181603	9597-GC-181604
18x18 in.	9597-GC-181802	9597-GC-181803	9597-GC-181804
20x16 in.	9597-GC-201602	9597-GC-201603	9597-GC-201604
20x18 in.	9597-GC-201802	9597-GC-201803	9597-GC-201804
22x18 in.	9597-GC-221802	9597-GC-221803	9597-GC-221804
24x18 in.	9597-GC-241802	9597-GC-241803	9597-GC-241804
Custom sizes available - any length, width or height.			

APEX Amputee Gel-Foam Cushion

RECOMMENDED USE: For amputee patients at risk of ulcers.

- Provides support and stability for post operative amputees
- MDF Wood base provides stability
- Available in left, right and bi-lateral



Size	Right	Left
16x16x3 in.	9597-AMP-161603-R	9597-AMP-161603-L
18x16x3 in.	9597-AMP-181603-R	9597-AMP-181603-L
18x18x3 in.	9597-AMP-181803-R	9597-AMP-181803-L
20x16x3 in.	9597-AMP-201603-R	9597-AMP-201603-L
20x18x3 in.	9597-AMP-201803-R	9597-AMP-201803-L
22x16x3 in.	9597-AMP-221603-R	9597-AMP-221603-L
22x18x3 in.	9597-AMP-221803-R	9597-AMP-221803-L
24x18x3 in.	9597-AMP-241803-R	9597-AMP-241803-L
26x18x3 in.	9597-AMP-261803-R	9597-AMP-261803-L
Custom sizes available - any length, width or height.		

APEX CORE™ Series

Meets CA 117 standards



The NYOrtho APEX CORE™ Gel-Foam Cushion combines a firm core foundation that eliminates hammocking with the superlative APEX™ innovative design that chambers the gel-pod closer to the top of the cushion, providing a cooler interface between the patient and wheelchair. A high-performance stretch cover improves moisture evaporation and enhances the APEX™ design features.

RECOMMENDED USE: Patients at risk for pressure ulcers.



APEX CORE Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers.



Size	Gel-Foam
16x16x3 in.	9599-G-161603
18x16x3 in.	9599-G-181603
18x18x3 in.	9599-G-181803
20x16x3 in.	9599-G-201603
20x18x3 in.	9599-G-201803
22x18x3 in.	9599-G-221803
24x18x3 in.	9599-G-241803
Custom sizes available - any length, width or height.	

APEX CORE Wedge Gel-Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring posture support, that are inclined to slide forward.



Size	Gel-Foam
16x16x4x2 in.	9599-GW-161642
18x16x4x2 in.	9599-GW-181642
18x18x4x2 in.	9599-GW-181842
20x16x4x2 in.	9599-GW-201642
20x18x4x2 in.	9599-GW-201842
22x18x4x2 in.	9599-GW-221842
24x18x4x2 in.	9599-GW-241842
Custom sizes available - any length, width or height.	

APEX CORE Pommel Gel-Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring hip positioning.



Size	Gel-Foam	Foam
16x16x3 in.	9599-GP-161603	9599-P-161603
18x16x3 in.	9599-GP-181603	9599-P-181603
18x18x3 in.	9599-GP-181803	9599-P-181803
20x16x3 in.	9599-GP-201603	9599-P-201603
20x18x3 in.	9599-GP-201803	9599-P-201803
22x18x3 in.	9599-GP-221803	9599-P-221803
24x18x3 in.	9599-GP-241803	9599-P-241803
Custom sizes available - any length, width or height.		

APEX CORE Wedge Pommel Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers requiring posture support, hip positioning and/or inclined to slide forward.



Size	Gel-Foam	Foam
16x16x4x2 in.	9599-GWP-161642	9599-WP-161642
18x16x4x2 in.	9599-GWP-181642	9599-WP-181642
18x18x4x2 in.	9599-GWP-181842	9599-WP-181842
20x16x4x2 in.	9599-GWP-201642	9599-WP-201642
20x18x4x2 in.	9599-GWP-201842	9599-WP-201842
22x18x4x2 in.	9599-GWP-221842	9599-WP-221842
24x18x4x2 in.	9599-GWP-241842	9599-WP-241842
Custom sizes available - any length, width or height.		

APEX CORE Coccyx Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure injuries in the coccyx area.



Size	Gel-Foam
16x16x3 in.	9599-GC-161603
18x16x3 in.	9599-GC-181603
18x18x3 in.	9599-GC-181803
20x16x3 in.	9599-GC-201603
20x18x3 in.	9599-GC-201803
22x18x3 in.	9599-GC-221803
24x18x3 in.	9599-GC-241803
Custom sizes available - any length, width or height.	

APEX CORE Coccyx Pommel Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers in the coccyx area, requiring hip positioning.



Size	Gel-Foam
16x16x3 in.	9599-GCP-161603
18x16x3 in.	9599-GCP-181603
18x18x3 in.	9599-GCP-181803
20x16x3 in.	9599-GCP-201603
20x18x3 in.	9599-GCP-201803
22x18x3 in.	9599-GCP-221803
24x18x3 in.	9599-GCP-241803
Custom sizes available - any length, width or height.	

Custom Cushions Upon Request

Any length, width or height.

APEX CORE Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers.



Size	Foam
16x16x3 in.	9599-161603
18x16x3 in.	9599-181603
18x18x3 in.	9599-181803
20x16x3 in.	9599-201603
20x18x3 in.	9599-201803
22x18x3 in.	9599-221803
24x18x3 in.	9599-241803
Custom sizes available - any length, width or height.	

APEX CORE Wedge Foam Cushion

RECOMMENDED USE: Patients at low risk for pressure ulcers, requiring posture support, that are inclined to slide forward.



Size	Foam
16x16x4x2 in.	9599-W-161642
18x16x4x2 in.	9599-W-181642
18x18x4x2 in.	9599-W-181842
20x16x4x2 in.	9599-W-201642
20x18x4x2 in.	9599-W-201842
22x18x4x2 in.	9599-W-221842
24x18x4x2 in.	9599-W-241842
Custom sizes available - any length, width or height.	

APEX CORE Anti-Thrust Gel-Foam Cushion Anti-Thrust Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers requiring anterior buildup at mid-line while preventing forward ischial migrations.



Size	Gel-Foam	Foam
16x16x4x2 in.	9599-ATG-161642	9599-ATH-161642
18x16x4x2 in.	9599-ATG-181642	9599-ATH-181642
18x18x4x2 in.	9599-ATG-181842	9599-ATH-181842
20x16x4x2 in.	9599-ATG-201642	9599-ATH-201642
20x18x4x2 in.	9599-ATG-201842	9599-ATH-201842
22x18x4x2 in.	9599-ATG-221842	9599-ATH-221842
24x18x4x2 in.	9599-ATG-241842	9599-ATH-241842
Custom sizes available - any length, width or height.		

APEX Anti-Thrust Foam Cushion Anti-Thrust Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers, requiring anterior buildup at mid-line while preventing forward ischial migrations.



Size	Gel-Foam	Foam
16x16x4x2 in.	9597-ATG-161642	9597-ATH-161642
18x16x4x2 in.	9597-ATG-181642	9597-ATH-181642
18x18x4x2 in.	9597-ATG-181842	9597-ATH-181842
20x16x4x2 in.	9597-ATG-201642	9597-ATH-201642
20x18x4x2 in.	9597-ATG-201842	9597-ATH-201842
22x18x4x2 in.	9597-ATG-221842	9597-ATH-221842
24x18x4x2 in.	9597-ATG-241842	9597-ATH-241842

Custom sizes available - any length, width or height.

Custom Cushions Upon Request

Any length, width or height.

APEX No Slide Cushion



Size	Gel-Foam
16x16x3 in.	9597-NSG-161603
18x16x3 in.	9597-NSG-181603
18x18x3 in.	9597-NSG-181803
20x16x3 in.	9597-NSG-201603
20x18x3 in.	9597-NSG-201803
22x18x3 in.	9597-NSG-221803
24x18x3 in.	9597-NSG-241803

Size	Foam
16x16x3 in.	9597-NSG-161603
18x16x3 in.	9597-NSG-181603
18x18x3 in.	9597-NSG-181803
20x16x3 in.	9597-NSG-201603
20x18x3 in.	9597-NSG-201803
22x18x3 in.	9597-NSG-221803
24x18x3 in.	9597-NSG-241803

APEX Visco Foam Cushion



Size	Visco-Foam
16x16x3 in.	9597-VSC-16163Q
18x16x3 in.	9597-VSC-18163Q
18x18x3 in.	9597-VSC-18183Q
20x16x3 in.	9597-VSC-20163Q
20x18x3 in.	9597-VSC-20183Q
22x18x3 in.	9597-VSC-22183Q
24x18x3 in.	9597-VSC-24183Q

APEX Waterfall Removable Pommel Foam Cushion



Size	Foam
16x16x3 in.	9596-RP-161603
18x16x3 in.	9596-RP-181603
18x18x3 in.	9596-RP-181803
20x16x3 in.	9596-RP-201603
20x18x3 in.	9596-RP-201803

WaveForm™ Series

Meets CA 117 standards



The **NYOrtho WaveForm Gel-Foam Cushion** is designed to complement the body's natural form and movement. As the name implies, the waveform-inspired design conforms to the body's shape more effectively, providing better support. A waterfall front curves to the natural shape of the leg trough while lateral and medial leg contours promote a slight leg abduction. A coccyx cutout reduces pressure in the affected area.

ABDUCTION

Medial thigh separator positions the legs, providing slight abduction for enhanced comfort and stability.

COCCYX CUTOUT

Coccyx cutout alleviates pressure on the tail-bone area, increasing comfort during prolonged periods of

WATERFALL FRONT

Waterfall front slopes downward, offering additional comfort for the legs and contributing to a more natural seated position.

DUAL GEL-POD

Internal gel-pod minimizes gel migration in the area of highest pressure and dissipates metabolic heat, continuously maintaining a cool interface between the body and cushion.

VISCO

Dual density foam layers provides a soft visco top entry layer promoting direct interface with cooling gel-pod and a higher density base layer for support and pressure distribution.

APEX Waveform Gel-Foam Cushion

RECOMMENDED USE: Patients at risk for pressure ulcers, requiring posture support and/or hip positioning.



SIZE	3 INCH
16x16x3 in.	9597-GLC-161603
18x16x3 in.	9597-GLC-181603
18x18x3 in.	9597-GLC-181803
20x16x3 in.	9597-GLC-201603
20x18x3 in.	9597-GLC-201803
22x18x3 in.	9597-GLC-221803
24x18x3 in.	9597-GLC-241803
Custom sizes available - any length, width or height.	

APEX Waveform Wedge Gel-Foam Cushion

Wedge design prevents forward sliding.

RECOMMENDED USE: Patients at risk for pressure injuries, requiring posture support, that are inclined to slide forward.



SIZE	
16x16x4x2 in.	9597-GLCW-161642
18x16x4x2 in.	9597-GLCW-181642
18x18x4x2 in.	9597-GLCW-181842
20x16x4x2 in.	9597-GLCW-201642
20x18x4x2 in.	9597-GLCW-201842
22x18x4x2 in.	9597-GLCW-221842
24x18x4x2 in.	9597-GLCW-241842
Custom sizes available - any length, width or height.	

APEX Waveform Anti-Thrust Gel-Foam Cushion

Anti-thrust design prevents forward ischial migration.

RECOMMENDED USE: Patients at risk for pressure ulcers requiring anterior buildup at mid-line to prevent forward ischial migrations.



SIZE	
16x16x4x2 in.	9597-CATG-161642
18x16x4x2 in.	9597-CATG-181642
18x18x4x2 in.	9597-CATG-181842
20x16x4x2 in.	9597-CATG-201642
20x18x4x2 in.	9597-CATG-201842
22x18x4x2 in.	9597-CATG-221842
24x18x4x2 in.	9597-CATG-241842
Custom sizes available - any length, width or height.	

TRiForm™ Ultra

Meets CA 117 standards



The **NYOrtho TRiForm Ultra Gel** is engineered with ultra-high density gel-infused memory foam to address specific pressure points and postural issues, incorporating an ischial step, coccyx scoop, and medial thigh separator to optimize pressure distribution, pelvic stability, and hip alignment. By addressing these biomechanical aspects, the cushion significantly reduces the risk of pressure ulcers and musculoskeletal complications. Crafted for long-term wheelchair users and those with particular medical needs, the cushion's multi-faceted design offers a scientifically-grounded approach to enhancing both immediate comfort and long-term musculoskeletal health.

RECOMMENDED USE: Patients at risk for pressure ulcers under the thighs, greater trochanters, ischial/tail-bone region and coccyx.



SIZE	
16x16x3x2 in.	9593E-161632
18x16x3x2 in.	9593E-181632
18x18x3x2 in.	9593E-181832
20x16x3x2 in.	9593E-201632
20x18x3x2 in.	9593E-201832
22x18x3x2 in.	9593E-221832
24x18x3x2 in.	9593E-241832
Custom sizes available - any length, width or height.	

MEDIAL THIGH SEPARATOR


PREVENTS EXCESSIVE INTERNAL ROTATION OF THE THIGHS.

COCCYX SCOOP

TARGETED RELIEF FOR TAIL BONE PRESSURE

ISCHIAL STEP

SUPPORTS ISCHIAL TUBEROSITIES, THE "SITTING BONE".



ANTI-HAMMOCKING DESIGN

PREVENTS SAGGING AND UNEVEN PRESSURE DISTRIBUTION

ULTRA-HIGH DENSITY GEL INFUSED MEMORY FOAM

CONTOURS SEAMLESSLY, PROVIDES CONSISTENT COOLING & UNIFORM FIRMNESS.

BENEFITS OF ULTRA-HIGH DENSITY GEL INFUSED MEMORY FOAM

Cooling: Ultra-high density gel-infused memory foam integrates gel throughout the foam structure, allowing for more even heat dispersion and consistent cooling across the entire surface, whereas a foam cushion with a separate gel-pack might have localized cooling only where the gel-pack is placed.

Consistent Firmness: Gel infusion can enhance the foam's responsiveness without drastically changing its firmness, ensuring a consistent feel throughout. A separate gel-pack could create inconsistencies in the cushion's firmness and feel.

Soft Contouring: Gel-infused memory foam still offers the contouring benefits of traditional memory foam, molding to the body's shape for personalized comfort.

Integrated Design: The infusion process blends the gel and foam seamlessly, eliminating the potential discomfort or unevenness that might arise from a separately placed gel-pack in a foam cushion.

NYOrtho Comfort Seats

Meets CA 117 standards

NYOrtho comfort seats provide protection of sensitive skin and comfortable seating. Comfort seats are secured in place with adjustable velcro straps or elastic and designed for use with positioning belts and many other accessories. Available in many different sizes, all comfort seats are machine washable and offer many different fabric options to suit your needs.

RECOMMENDED USE: Patients at risk for pressure ulcers.



METRIS™ Comfort Seat MODEL 9522GC-GEL

Actual Dimensions: 72"L x 19"W x 3"H



Easy-To-Attach
Quick-Release Straps

Visco Elastic Memory Foam

Visco memory foam supports and contours upper body shape for maximum comfort.

Rotational Cover Fluid-Proof

SURE-CHEK™ thermoplastic fabric composites are strong, durable, antimicrobial, anti-static, flame resistant, stain, odor, and fluid proof.

Rotational Cover Fluid-Resistant

Super soft **PEACHSKIN** fabric is strong, durable, antimicrobial, and fluid resistant - improves moisture evaporation. CA 117-2013

Independent-Cell Foam Body

High-performance independent-cell foam cavity provides greater air-flow and circulation to the legs.



ULTRA-LIFT™ Comfort Seat MODEL 9523GC-GEL

Actual Dimensions: 65"L x 19"W x 3"H



Easy-To-Attach
Quick-Release Straps

Visco Elastic Memory Foam

Visco memory foam supports and contours upper body shape for maximum comfort.

Foam / Gel Seating Surface

High-density foam surrounds an internal dual-chamber gel pod, providing a cool and comfortable interface that evenly maintains gel in the area of the highest pressure, preventing skin breakdown.

Fluid-Proof Cover

SURE-CHEK™ thermoplastic fabric composites are strong, durable, antimicrobial, anti-static, flame resistant, stain, odor, and fluid proof.

Independent-Cell Heel Suspension Cushion

Wedge shaped foam cushion elevates heel providing zero-g pressure.



Geri-Chair Gel-Foam Comfort Seat

MODEL 9525-2G-021970

Actual Dimensions: 72"L x 19"W x 3"H



Easy-To-Attach
Quick-Release Straps

Foam / 4Gel Seating Surface

High-density foam surrounds an internal quad-chamber gel pod, providing a cool and comfortable interface that evenly maintains gel in the area of the highest pressure, preventing skin breakdown.

Breathable, Stretch Cover

High performance pressure reducing fabric promotes healing via active evaporation of moisture and cooling at patient pressure points.



Geri-Chair Comfort Seat

RECOMMENDED USE: Patients requiring long-term seating comfort and/or protection of sensitive skin.



9522GC



9522GC-S

Model	Style	Color	Size
9519GC	Antimicrobial / Water-Resistant	Light Blue	Full 71x18 in.
9519GC-S	Antimicrobial / Water-Resistant	Light Blue	Short 47x18 in.
9520GC	Water-Resistant	Navy	Full 71x18 in.
9520GC-S	Water-Resistant	Navy	Short 47x18 in.
9521GC	Soft, Breathable Cotton	*Fun	Full 71x18 in.
9521GC-S	Soft, Breathable Cotton	*Fun	Short 47x18 in.
9522GC	Sure-Chek/WR Rotational Cover	Gray / Light Blue	Full 71x18 in.
9522GC-S	Sure-Chek/WR Rotational Cover	Gray / Light Blue	Short 47x18 in.
Custom sizes available - any length or width.		*Style may change periodically.	

Multiple Fabric Choices:

- Super-soft peachskin fabric: anti-microbial / water-resistant
- Low shear taslon fabric: water-resistant
- Designed with fun cotton fabric – a soft, breathable and durable fabric with excellent wickability.
- Rotational Cover: **SURE-CHEK®** thermoplastic fabric composites are strong, durable, antimicrobial, anti-static, fluid proof, and flame, stain, and odor resistant. Other side: Super soft **PEACHSKIN** fabric is strong, durable, antimicrobial, and fluid resistant – improves moisture evaporation.

Wheelchair Comfort Seat

RECOMMENDED USE: Patients requiring long-term seating comfort and/or protection of sensitive skin.



9521

Model	Style	Color	Size
9519	Antimicrobial / Water-Resistant	Light Blue	Short 36x17 in.
9519L	Antimicrobial / Water-Resistant	Light Blue	Reg 47x17 in.
9519XL	Antimicrobial / Water-Resistant	Light Blue	Full 54x17 in.
9520	Water-Resistant	Navy	Short 36x17 in.
9520L	Water-Resistant	Navy	Reg 47x17 in.
9520XL	Water-Resistant	Navy	Full 54x17 in.
9521	Soft, Breathable Cotton	*Fun	Short 36x17 in.
9521L	Soft, Breathable Cotton	*Fun	Reg 47x17 in.
9521XL	Soft, Breathable Cotton	*Fun	Full 54x17 in.
9522	Sure-Chek/WR Rotational Cover	Gray / Light Blue	Short 36x17 in.
9522L	Sure-Chek/WR Rotational Cover	Gray / Light Blue	Reg 47x17 in.
9522XL	Sure-Chek/WR Rotational Cover	Gray / Light Blue	Full 54x17 in.
Custom sizes available - any length or width.		*Style may change periodically.	



9519L



9520XL

- All extensions are removable with a convenient zipper feature on the 9519, 9520 and 9521 styles only.

WHEELCHAIR ACCESSORIES

NYOrtho Wheelchair & Geri-chair accessories provide for resident protection and comfortable seating.

Wheelchair Lap Cushions

Meets CA 117 standards

NYOrtho Wheelchair Lap Cushions help assist in upper body alignment, providing upper torso support and an elevated resting spot for hands. Easy to clean vinyl cover.

RECOMMENDED USE: Patients needing upper torso support or an elevated resting spot for hands.



Self-Releasing



Full Arm



Desk Arm



- Ideal for residents with minimal hand/arm movement, releases automatically as the resident rises.
- Self-releasing model designed with a simple-to-use hook and loop closure for easy self-release. #9526
- Full arm design designed with front notches to accommodate a full arm wheelchair. #9527
- Desk arm design has side notches to fit desk arm wheelchairs. #9528

Model	Style	Dimensions	Fits W/C
9526	Self-Releasing	4x9x20 in.	18-22 in.
9527-1620-3	Full Arm	3x10x27 in.	16-20 in.
9527-2224-3	Full Arm	3x10x31 in.	22-24 in.
9528-18-3	Desk Arm	3x8x27 in.	18-22 in.
Custom sizes available - any length, width or height.			

Smoker's Apron

RECOMMENDED USE: Designed to assist facilities in providing safeguard for patients from smoking materials and flame.



Ideal for patients requiring a protective covering to shield against a burning match or lit cigarette while accommodating patient/resident rights.

- Certified flame-resistant fabric meets **ASTM F2302**.
NYOrtho Smoker's Apron meets compliance requirements for use in institutional facilities, such as nursing homes, assisted living, and similar facilities.
- Flame-resistant apron protects against smoldering fires and cannot be ignited by a cigarette, cigar or burning match.
- Side straps prevent cigarettes and ashes from falling between resident and the wheelchair.
- Wipe clean, lightweight and foldable.
- Universal sizing.

Model	Style	Color	Dimensions	Fits
9530	Standard	White	30x34 in.	W/C
9531	Standard	White	30x44 in.	Geri-Chair
9532	Deluxe - Exceeds ASTM F2302	Gray	30x34 in.	W/C

Wheelchair Lateral Supports

Meets CA 117 standards

NYOrtho Wheelchair Lateral Supports provide cushioned support while correcting lateral leaning, keeping the upper torso aligned and correcting sitting posture.



Wheelchair Deluxe Lateral Arm Support

RECOMMENDED USE: For patients who tend to lean laterally or who need arm elevation.



The **NYOrtho Deluxe Lateral Arm Support** helps correct lateral leaning while providing a comfortable foam filled base to rest the arm. Features a high-performance wipe-clean, stretch cover that improves moisture evaporation.

- Contoured design comfortably cradles the arm while providing arm elevation to reduce lateral leaning.
- Double straps maintain arm position to reduce lateral leaning.
- Wipe-clean, high-performance pressure reducing fabric. Antimicrobial, water-resistant and mildew-resistant.
- Can be used on either wheelchair arm.



Model	Dim.	Fits	UOM
9533	9x16 in.	Most Armrests	Each

Wheelchair Lateral Stabilizer Armrest Bolster

RECOMMENDED USE: Patients who lean laterally.



The **NYOrtho Lateral Stabilizer Armrest Bolster** provides cushioned support while helping control upper torso alignment.

- Comfortable arm support promotes better posture and stabilizes resident.
- Available in a cool and comfortable, high-performance stretch fabric.
- Available in full arm or desk arm length.
- Flat top cushion keeps arm in place and greatly improves patient stability.
- Round top cushion provides body and arm cushioning and controls upper torso alignment.
- Easy-on, easy-off slit opening easily fits over armrest and secures with hook and loop closure.
- Wipe-clean surface.

Model	Dim.	Style	UOM
9534D-FD-0611	6x11 in.	Flat Desk Arm	Each
9534D-FF-0616	6x16 in.	Flat Full Arm	Each
9534D-RD-0611	6x11 in.	Round Desk Arm	Each
9534D-RF-0616	6x16 in.	Round Full Arm	Each



Wheelchair Lateral Support L-Shape

RECOMMENDED USE: Patients who require bilateral upper torso support.



The **NYOrtho Wheelchair Lateral Support L-Shape** provides petite residents with bilateral torso support and promote proper sitting posture in a standard wheelchair or geri-chair. Made with MDF rigid insert, covered with padded foam and high-performance wipe-clean material.

- Improves comfort, filling the gap between resident and armrest.
- MDF rigid insert promotes better posture and stabilizes patient position.
- Cool and comfortable, high-performance stretch fabric covers an inner foam layer.
- Easy to apply reversible elastic bands keep support in place.
- Sold as an each or as a pair.
- Wipe-clean surface.

Model	Dim.	UOM
9535-16131H	13x16 in.	Each
9535-16131H-S	13x16 in.	Pair

Wheelchair BalancePro Lateral Support

RECOMMENDED USE: For patients who tend to lean laterally or who need arm elevation.



The NYOrtho BalancePro Lateral Support addresses various physical conditions that compromise normal muscular strength, control, or function, when the body's ability to maintain upright posture or evenly distribute weight is compromised, leading to a tendency to lean to one side. This can be due to muscle weakness, muscle loss, uneven muscle strength, paralysis, skeletal deformities, improper wheelchair fit, or poor postural habits.

- Raised, trough shaped arm cushion maintains arm position and upright position.
- Wipe-clean, high-performance pressure reducing fabric. Antimicrobial, water-resistant and mildew-resistant.
- Can be used on either wheelchair arm.
- Wipe-clean surface.
- CA 117-2013.

Model	Dim.	UOM
9537	Universal	Each

Lateral Support Assembly 3-Part System

RECOMMENDED USE: For patients at risk of lateral leaning when seated in a wheelchair or geri-chair.



The **NYOrtho Lateral Support Assembly 3-Part System** is designed to improve patient's position and comfort when seated in a wheelchair or geri-chair. Twin lateral foam wings, adjustable in any direction, are covered with a high-performance stretch cover that improves moisture evaporation.



- Lateral wings provide support, increase comfort and help correct lateral leaning when seated.
- Easy-to-install 3-part system fits over existing sling back.
- Lateral wings attach to back with hook and loop attachment, allowing for quick and easy adjustments in any direction.
- High-performance pressure reducing fabric promotes healing via active evaporation of moisture and cooling at patient pressure points. Antimicrobial, water-resistant and mildew-resistant.
- Fits wheelchairs and geri-chairs.
- Wheelchair models available in 4" or 6" depth pads.

Model	Fits	Style	UOM
9543-1622	W/C 16-22 in.	Wheelchair - 4" Pad	Set
9543-1622-6	W/C 16-22 in.	Wheelchair - 6" Pad	Set
9543G-1622	G/C 16-22 in.	Geri-Chair - 4" Pad	Set

CurveFit Wheelchair Support

RECOMMENDED USE: Patients with spinal conditions that alter the spine's natural curvature.

The **NYOrtho CurveFit Wheelchair Support** offers a customizable backrest for wheelchairs, adjusting its curvature to align with the patient's unique spinal shape, providing tailored support. This design not only reduces pain and discomfort by ensuring even weight distribution but also encourages extended sitting, benefiting those attending events, traveling, or working. Moreover, it can aid respiratory function in severe kyphosis cases and foster increased social interaction due to enhanced comfort. Being an attachable accessory, it's a cost-effective solution, eliminating the need for an entirely new custom wheelchair.



Model	Fits	UOM
8755	16-22 in.	Each

- **Ergonomic Contoured Shape:** Cushion angles away from the frame, eliminating any contact with wheelchair frame.
- **Effortless Tailored Support** Incorporating hook and loop straps to customize its curvature to align with the patient's unique spinal shape.
- **High-Performance Fabric:** pressure reducing fabric promotes healing via active evaporation of moisture and cooling at patient pressure points. Antimicrobial, water-resistant and mildew-resistant.
- **Built to Last:** Durable construction guarantees longevity and unwavering protection over time.
- **Safety:** Meets CA 117-2013 standards.



Conditions:

Various spinal conditions alter the spine's natural curvature, leading to discomfort and abnormal posture. **Kyphosis**, characterized by an excessive outward curve, and **scoliosis**, a sideways curvature, can either be congenital or develop later in life. **Lordosis** involves an exaggerated inward lumbar curve. **Ankylosing spondylitis** is an arthritis type that results in vertebrae fusion, leading to a forward-stooped posture. Consult a healthcare professional before making mobility and support equipment decisions.

Foot & Calf Pad

RECOMMENDED USE: Patients needing lower extremity support and/or cushioned leg support while preventing feet from slipping behind footrest.

The **NYOrtho Foot & Calf Pad** is tailored to provide exceptional comfort while safeguarding delicate skin from potential harm in the lower extremities. This accessory envelops the wheelchair's footrests and leg rest pads, acting as a protective barrier against skin tears, bruises, and related injuries. Attaching it is a breeze thanks to its quick-release buckles and reliable hook and loop fasteners on both sides, ensuring a snug and durable fit.

Shielding the Essential Danger Zones: Comprehensive protection from both footrest and leg rest panels.

- **Ergonomic Contoured Shape:** Cushion angles away from the frame, eliminating any contact with wheelchair frame.
- **Effortless Tailored Support** Incorporating hook and loop straps to customize its curvature to align with the patient's unique spinal shape.
- **High-Performance Fabric:** pressure reducing fabric promotes healing via active evaporation of moisture and cooling at patient pressure points. Antimicrobial, water-resistant and mildew-resistant.
- **Built to Last:** Durable construction guarantees longevity and unwavering protection over time.
- **Safety:** Meets CA 117-2013 standards.

Model	Fits	UOM
8750	16-20 in.	Each



Wheelchair Urinary Drain Bag Holder

RECOMMENDED USE: Patients with urine drainage bags.



Urinary Drain Bag Holders provide a convenient and modest covering of urinary drainage bags for when residents are catheterized.

- Preserves resident dignity by hanging discreetly under wheelchair, geri-chair or bed.
- Secured in place with adjustable straps.
- Universally sized, holds all drainage-bags.
- Available in wipe-clean vinyl w/hook and loop closure or machine washable canvas w/sliding buckle.
- Drain bag not included.

Model	Style	Fits W/C	UOM
9549	Canvas	Universal	Each
9550	Vinyl	Universal	Each



Wheelchair Armrest Covers

RECOMMENDED USE: Patients requiring protection of sensitive skin and/or improved seating comfort.



NYOrtho Wheelchair Armrest Covers are designed to protect sensitive arms and elbows against contact with hard armrests and metal wheelchair frames. Covers are made with high-performance stretch fabric which remains cool and soothing to the skin.

- Cool and comfortable, high-performance stretch fabric covers an inner foam layer.
- Provides a comfortable base, protecting sensitive skin while improving seating comfort.
- Easy to attach hook and loop closure allows cover to be reversed based on your color choice.
- Sold as a pair.
- Wipe-clean surface.

Model	Dim.	UOM
9536-11	11x11 in.	Pair
9536-16	11x16 in.	Pair
9536-18	11x18 in.	Pair
9536-22	11x22 in.	Pair

Reversible



Abduction Wedge Knee Separator

RECOMMENDED USE: Patients with abduction contractures, recovering from hip replacement surgery and/or at risk of skin breakdown.



NYOrtho Abduction Wedge Knee Separator is most commonly used to keep the legs stable after hip replacement surgery or for abduction contractures.

- Contoured sides maintain 4" knee separation.
- Available in two fabric options: light gray Sure-Chek fabric – Antimicrobial, fluid-proof, hypoallergenic, wipe-clean surface, or royal blue, breathable, stretch fabric – Antimicrobial.
- Hook receptive medical foam strap comfortably maintains leg position, excess strap can be trimmed.
- Strap is 46" / Foam 4".
- For wheelchair or bed.
- Meets CA 117.

Model	Style	UOM
9820-96	Sure-Chek Cover	Each
9820-97	Breathable Stretch	Each



Breathable Stretch Cover

Urine Bottle Holder

RECOMMENDED USE: Provides a secure holder for patients requiring urine bottles when in a wheelchair, geri-chair or bed.



The **NYOrtho Urine Bottle Holder** provides a convenient, discreet and sanitary pouch for holding a male urinal. Male patients using urine bottles will find this holder provides a discreet way to hold a urine bottle and keep it easily accessible.

- Easy-on easy-off quick-release buckle attaches to wheelchairs and bed frames.
- Gripper webbing keeps holder in place and is easy to adjust.
- Stiff stay maintains shape for easy removal and replacement of holder.
- Durable and washable nylon.
- For use with 1000ml male urine bottle.

Model	Fits	UOM
9544	1000ml male urine bottle	Each

Urine bottle not included



Wheelchair Oxygen Tank Holder

RECOMMENDED USE: Provides a secure holder for patients requiring oxygen tanks when in a wheelchair.



For patients on the go that require oxygen, the **NYOrtho Wheelchair Oxygen Cylinder Holder** safely and comfortably stores your patient's oxygen tank off center, behind the wheelchair back, to avoid interference with the wheelchair user's head or back. Adjustable straps accommodate all wheelchair widths.

- Economical alternative to steel tank holders.
- Adjustable straps position the tank off center to avoid patient's head or back.
- Made of easy-to-clean, heavy-duty, waterproof nylon.
- Attaches with easy-to-adjust quick-release buckles.
- Dimensions: 27 3/4" L x 7 1/2" W.
- Fits Size D (M15) and E (M24) cylinders. Urine bottle not included.

Model	Tank Size	Fits W/C	UOM
9546	D&E	16-22 in.	Each
9546-2430	D&E	24-30 in.	Each



Wheelchair Foot & Leg Rest Pads

Meets CA 117 standards



The NYOrtho Footrest Extender Leg Rest Pad extends the base of support and elevates the legs one, two or three inches, improving posture, promoting comfort and alleviating foot drop symptoms. A cushioned support behind the leg helps prevent the feet from slipping behind the footrest and getting caught under the wheelchair. Covered in durable, wipe-clean vinyl.

Wheelchair Leg Rest Pad



- Cushioned Leg Rest Pad protects lower leg against contact with wheelchair and prevents feet from slipping behind footrests.
- Secures easily with quick-release straps.



Model	Dim.	Fits W/C	UOM
9547	16x9 in.	16-20 in.	Each
9547-22	22x9 in.	22-24 in.	Each
Custom sizes available - any length, width or height.			

Wheelchair Footrest Extender Leg Rest Pad



- Modular footrest extender raises footrest 1, 2 or 3 inches minimizing foot drop.
- Cushioned Leg Rest Pad protects lower leg against making uncomfortable contact with wheelchair and prevents feet from slipping behind footrests.
- Secures easily with quick-release straps.

Model	Style	Dim.	Fits W/C	UOM
9548-1620-1	1" Extender	15.75x8.25 in.	16-20 in.	Each
9548-1620-2	2" Extender	15.75x8.25 in.	16-20 in.	Each
9548-1620-3	3" Extender	15.75x8.25 in.	16-20 in.	Each
9548-2224-1	1" Extender	25.75x8.25 in.	22-24 in.	Each
9548-2224-2	2" Extender	25.75x8.25 in.	22-24 in.	Each
9548-2224-3	3" Extender	25.75x8.25 in.	22-24 in.	Each
Custom sizes available - any length, width or height.				



Wheelchair Footrest Bag



- Provides storage of personal items and is a perfect stowaway for the Footrest Extender and Leg Rest Pad.
- Easily hangs off wheelchair.
- Covered in durable, wipe-clean vinyl.



Model	Dim.	Fits W/C	UOM
9548B-1620	16x23 in.	16-20 in.	Each
9548B-2226	22x23 in.	22-26 in.	Each
Custom sizes available - any length or height.			

CushionGuard Wheelchair Footrest Sleeve

RECOMMENDED USE: For patients who need to maintain their feet in a stable position and require a padded footrest.



The **NYOrtho CushionGuard Wheelchair Footrest Sleeve** offers an essential combination of comfort, safety, and hygiene for wheelchair users. Its foam-filled, non-slip design ensures user comfort and stability, while the antimicrobial fabric promotes cleanliness and health. Easy to apply, clean and durable, this footrest cover also allows personal customization, enhancing the wheelchair experience.



Model	Dim.	UOM
9542	6.75 x 6.75 in.	Each

Sheepskin Wheelchair Pads

RECOMMENDED USE: Patients requiring protection of sensitive skin and/or improved seating comfort.



Soft, plush synthetic sheepskin transforms any wheelchair, adding pliable comfort while protecting skin from vinyl surface.

- High-pile fur promotes air circulation beneath the resident and provides comfortable cushioning.
- Low shear, non-allergenic synthetic sheepskin reduces risk of pressure sores.
- Provides a comfortable base preventing skin from sticking to wheelchair vinyl seats and backs.
- Easy to attach hook and loop closure, elastic or straps – depends on the model.
- Machine washable and meets CA 117 Standards.
- Berry Compliant Fabric.
- Sold as a set or individually.



Model	Style	Dim.	UOM
9552	Complete Set		Set
9553	Armrest Pads - Hook & Loop	8.5x15 in.	Pair
9554	Leg Pad - Adjustable Strap	12x23 in.	Each
9555	Footrest Pads - Elastic	7x7 in.	Pair
9556	Seat & Backrest Pads - Ties	18x16 in.	Set

Custom sizes available - any length or height.



BED SAFETY & ACCESSORIES

NYOrtho bed safety products and accessories provide safety and comfort for your patients.

Perimeter Mattress Cover

RECOMMENDED USE: Patients only needing a gentle reminder of the location of the mattress edges, and help provide a less restrictive environment than side rails.



The **NYOrtho Perimeter Mattress Cover** provides a less restrictive defined perimeter to help prevent falls. Complies with latest restraint reduction legislation to reduce the use of side rails, avoiding entrapment. An open mid section accommodates easy patient egress and ingress.

- Quick-release buckles easily secure the covers to the bed.
- Medium density foam bolsters are removable from zippered compartments to facilitate easy laundering of cover.
- Easy to remove foam foot bolster.
- Convenient sheet attachment clips keeps cover sheet in place.
- Nylon top cover is fluid resistant, vapor permeable and anti-microbial.
- Machine washable.
- Meets CA 117 standard.



Model	Dim.	UOM
9925-36	36 in.	Each
9925-39	39 in.	Each
9925-42	42 in.	Each
Custom sizes available - any length, width or height.		

Padded Wall Bumper

RECOMMENDED USE: For protection of walls and furniture when utilizing a hospital bed.



The **NYOrtho Padded Wall Bumper** protects wall and headboards from the inevitable damage that may occur in the course of use in hospitals or nursing homes. Easy to install on the headboard, the bumper is easily removed to wipe clean.

- Protects walls from damage by bed/headboard.
- Nylon top cover is fluid resistant, vapor permeable and anti-microbial.
- Easily attaches to headboard with adhesive backed hook and loop.
- Easy to remove.



Model	Dim.	UOM
9568-3H0638	3.5x6x38 in.	Each

Repositioning Slide

RECOMMENDED USE: Patients requiring assistance to reposition and at risk of skin breakdown/ friction burns.



The **NYOrtho Re-positioning Slide** sheet is used for re-positioning, turning and moving residents up in bed. The tubular design and multiple handles facilitate easy re-positioning, enabling two assistants to slide the patient in the direction of the headboard. A low-friction nylon fabric makes re-positioning easy and effortless, preventing caregiver neck and back injury and protects resident's skin integrity and allows for safer re-positioning by significantly reducing friction and shearing.

- Made from durable 400 denier bullet nylon.
- Machine washable.

Model	Dim.	UOM
9580-26H75	26.5x75 in.	Each
9580-3475	34x75 in.	Each

Bed Rail Pads

Meets CA 117 standards



NYOrtho Bed Rail Pads provide protection from bed rails. For extra resident protection, bed rail pads extend between and below the mattress and bed rails. Removable with easy to apply with hook and loop attachment. Rails may be raised or lowered without removing bed rail pads. Available in half and full rail sizes. Vinyl cover wipes clean for easy disinfecting.

RECOMMENDED USE: Patients at risk for bed entrapment; patients prone to unassisted bed exit.

Sheepskin Bed Rail Pads



Soft synthetic **Sheepskin Bed Rail Pads** provide comfort and protection to sensitive skin. Helps minimize accidental patient injuries that result from striking bedside rails.

- Provides protective padding to help minimize the risk of patient injury.
- Low shear, non-allergenic, synthetic sheepskin covers side and top of bed rail.
- Rails may be raised or lowered without removing bed rail pad.
- Pad is secured to rail with easy hook and loop closures and ties.
- Accommodates both half and full-size side rails.
- Machine washable.

Model	Dim.	UOM
9557-24	18x24 in.	Each
9557-30	18x30 in.	Each
9557-34	18x34 in.	Each
9557-42	18x42 in.	Each
9557-48	18x48 in.	Each
9557-56	18x56 in.	Each
9557-60	18x60 in.	Each
9557-72	18x72 in.	Each

Custom sizes available - any length or height.

Cushioned Top Vinyl Bed Rail Pads



- Cushioned top of rail for added protection.

Model	Dim.	UOM
9565-011537-C	1x15x37 in.	Pair
Custom sizes available - any length or height.		

Split Rail Vinyl Bed Rail Pads



Model	Dim.	UOM
9565-011124	1x11x24 in.	Pair
9565-011134	1x11x34 in.	Pair
Custom sizes available - any length or height.		

Vinyl Bed Rail Pads



Model	Dim.	UOM
9565-011169	1x11x69 in.	Pair
Custom sizes available - any length or height.		



- Stretch hook & loop closure fasten pads securely to bed rails on "S" models.

Model	Dim.	UOM
9565-010928-S	1x9x28 in.	Set of 4
Custom sizes available - any length or height.		

Bed Rail Wedge

RECOMMENDED USE: Patients at risk for bed entrapment and injury via side rails; patients prone to unassisted bed exit.



NYOrtho Bed Rail Wedges help protect patients from injury by closing the gaps around the mattress perimeter and by cushioning the side rails. They also help to deter the patient's unassisted bed exit.

- Features padding to help prevent patient injury.
- Hook & loop closure fasten pads securely to bed rails.
- Nylon top cover is fluid resistant, vapor permeable and antimicrobial.
- Available in half and full rail sizes.
- Meets CA 117-2013 standard.

Model	Dim.	UOM
9565W-35	35 in.	Pair
9565W-70	70 in.	Each
Custom sizes available - any length or height.		



Headboard Pad Wall Protector

RECOMMENDED USE: To close the gaps at the headboard of a bed; cushions head to help prevent injury.



The **NYOrtho Headboard Pad Wall Protector** is designed to prevent injury by cushioning the headboard and by closing the gap between the mattress and the headboard. Pad also functions as a wall protector.

- Features a padded bumper that covers the headboard to help prevent injury by closing the gaps in the bed.
- Also acts as a wall protector preventing damage from the bed continuously hitting the wall.
- Pad easily attaches to the headboard with straps and hook and loop closure.
- Nylon top cover is fluid resistant, vapor permeable and ant-microbial.

Model	Dim.	UOM
9566H-022236	2x22x36 in.	Each
Custom sizes available - any length or height.		



Gap Shield

RECOMMENDED USE: Patients at risk for bed entrapment.



NYOrtho Bed Rail Wedges help protect patients from injury by closing the gaps around the mattress perimeter and by cushioning the side rails. They also help to deter the patient's unassisted bed exit.

- Features padding to help prevent patient injury.
- Hook & loop closure fasten pads securely to bed rails.
- Nylon top cover is fluid resistant, vapor permeable and antimicrobial.
- Available in half and full rail sizes.
- Meets CA 117-2013 standard.

Model	Dim.	UOM
9567-35	35 in.	Pair
9567-42	42 in.	Each
9567-70	70 in.	Each
Custom sizes available - any length or height.		

Roll Control Bolsters

RECOMMENDED USE: Patients prone to exiting bed unassisted and/or falling out of bed.



A safe alternative to side rails, **NYOrtho Roll-Control Bolsters** help prevent falls from rolling out of bed. The adjustable bolster attaches to the bed with quick-release buckles and may be repositioned to accommodate individual patient needs. While providing sufficient space for patient comfort, the bolster is not considered a restraint for most patients. Roll control bolsters can be added to help reduce potential bed rail entrapment zones.

- Repositionable bolster attaches with quick-release buckle for individual positioning and is easily removed for patient transfer.
- Available in single unit.
- Covered in durable wipe clean vinyl.
- Meets CA 117 standard.

Model	Dim.	UOM
9912	7x8x34 in.	Set
Custom sizes available - any length or height.		

Mattress Extender

RECOMMENDED USE: Patients at risk for bed entrapment.



NYOrtho Mattress Extenders help eliminate gaps at the head or foot of the bed, providing a more comfortable and safe sleeping experience.

- Nylon top cover is fluid resistant, vapor permeable and anti-microbial.
- Medium density foam.
- Meets CA 117 standard.

Model	Dim.	UOM
9569-030636	3x6x36 in.	Each
9569-030642	3x6x42 in.	Each
9569-030648	3x6x48 in.	Each
9569-040636	4x6x36 in.	Each
9569-040642	4x6x42 in.	Each
9569-040648	4x6x48 in.	Each
9569-050636	5x6x36 in.	Each
Custom sizes available - any length or height.		

Lateral Position Knee Abductor

RECOMMENDED USE: For side-lying patients with leg contractures that require knee abduction.

The NYOrtho Lateral Position Knee Abductor is an expertly designed supportive device that aids in the management of leg contractures and promotes knee abduction. Specifically designed for use in a side-lying position on a bed, it provides optimal support to maintain the correct positioning of the lower limbs during rest or sleep. This device works by maintaining the proper alignment and separation of the knees and legs, which is essential in preventing and treating leg contractures. By promoting knee abduction, it not only ensures correct posture but also reduces the risk of pressure sores and aids in post-surgical recovery.



- **Prevention of Contractures:** Knee abductors can help prevent leg contractures by keeping the knees and legs properly aligned and separated.
- **Post-Surgical Recovery:** After knee surgery, especially after total knee replacement or procedures involving the lower limbs, a knee abductor can help maintain proper alignment and prevent the legs from touching each other, which can help in the healing process and prevent complications. Available in single unit.
- **Pressure Sore Prevention:** By keeping the knees slightly apart, a knee abductor can help reduce the pressure and friction that occurs when the knees rub together, which can lead to pressure sores, skin irritation, or exacerbate existing sores.
- **Pain and Discomfort Relief:** For individuals with certain conditions like arthritis or hip dysplasia, a knee abductor can help reduce discomfort and pain by promoting a more natural and comfortable alignment while resting or sleeping.
- **Improved Posture and Circulation:** By promoting proper alignment of the lower body, knee abductors can also improve overall body posture and enhance circulation in the lower limbs.

Model	Fits	UOM
9823	Universal	Ea.

Endurance Pillows



NYOrtho Endurance Pillows are engineered to be comfortable and environmentally friendly as per the Safe Guard Barrier Systems® certification. All pillows are filled with a hypo-allergenic polyester fill and meet NFPA 702 fire-rating standards.

- Nylon top cover is fluid resistant, vapor permeable and anti-microbial.
- Medium density foam.
- Meets CA 117 standard.



Model	Dim.	UOM
9620-1824	18x24 in.	6/Case
9620-2026	20x26 in.	6/Case
9620-2030	20x30 in.	6/Case

Sheepskin Bed Pads

RECOMMENDED USE: Patients at low risk of skin breakdown



NYOrtho Sheepskin Bed Pads help make any bed cozy and comfortable with soft, plush, synthetic sheepskin materials. Used in many products, synthetic sheepskin is soft to the touch and feels wonderful against the skin.

- Low shear, non-allergenic synthetic sheepskin reduces risk of pressure sores.
- High-pile fur promotes air circulation beneath the resident and provides comfortable cushioning.
- Available in many sizes to suit your needs.
- Machine washable.
- Meets CA 117 standards.

Model	Dim.	UOM
9560-2430	24x30 in.	Each
9560-3040	30x40 in.	Each
9560-3060	30x60 in.	Each
9560-4880	48x80 in.	Each
9560-30360	30x360 in.	Roll
Custom sizes available - any length or height.		

Wedges

Meets CA 117 standards



Positioning Bed Wedges are utilized to reposition bed-ridden patients, helping eliminate pressure points over the bony prominence's to help reduce the occurrence of pressure ulcers. Elevating Bed Wedges are utilized to elevate the lower or upper body to help improve circulation, relieve pressure on sensitive areas, reduce swelling, as well as to improve breathing.

RECOMMENDED USE: Patients requiring additional positioning and support.

Positioning Bed Wedges



- Available in many different angles of incline that help to relieve pressure on hipbone and coccyx.
- Breathable, low shear, stretch cover with antimicrobial protection reduces friction.
- Slip resistant bottom prevents sliding.

Model	Incline	Dim.	UOM
9910-161105	25°	5x11x16 in.	Pair
9910-161207	30°	7x12x16 in.	Pair
9910-241207	30°	7x12x24 in.	Pair
9910-341207	30°	7x12x34 in.	Pair
9910-170808	45°	8x8x17 in.	Pair
9910-281082	40°	2x8x10x28 in.	Each
Custom sizes available - any length, width or height.			

Sure-Chek Positioning Bed Wedges

AM

CSA

HA

WP

TC

Sure-Chek
Sometimes you need the best cover available.

MODEL 9910V--XXXXXX



• FLAME RESISTANT

All Sure-Chek healthcare fabrics are flame resistant and self-extinguish when a source flame is removed.

• THERMAL CONDUCTIVITY

Sure-Chek healthcare fabrics all perform within a desirable range for patient comfort without sacrificing important technical performance benefits.

• HYPOALLERGENIC

Sure-Chek healthcare fabrics have been engineered to be hypoallergenic and comfortable for the user.

• ANTIMICROBIAL FABRIC PROTECTION

The antimicrobial agent helps protect the fabric from microbial attack and reduces undesirable organic odors.

• FLUID PROOF & MOISTURE VAPOR TRANSMISSION

All Sure-Chek healthcare fabrics are fluid and moisture proof, reducing the risk of contamination of the mattress from liquids or bodily fluids.

Elevating Bed Wedges



Model	Dim.	UOM
9911-242407	7x24x24 in.	Each
9910-242410	10x24x24 in.	Each
9911-242412	12x24x24 in.	Each
Custom sizes available - any length, width or height.		

- Provides elevation to improve position.
- Removable water-resistant cover is machine washable.

Abduction Wedge

RECOMMENDED USE: Designed for use following hip surgery when immobilization or post-operative positioning is required.



The **NYOrtho Abduction Wedge** maintains and immobilizes following hip surgery or other post operative positioning.

- High-performance stretch cover that improves moisture evaporation.
- Wide polypro straps keeps positioner in place and allows for easy adjustment.
- Medium density foam.
- Meets CA 117 standard.

Model	Dim.	UOM
9822-18	6x12x18 in.	Each
9822-22	6x15x22 in.	Each
9822-25	6x18x25 in.	Each

Hip Abductor and Knee Extender Combo

RECOMMENDED USE: For individuals who may be bedridden, post-surgery patients, or those suffering from certain types of lower body injuries or ailments.



- **Hip Abduction:** The hip abductor component is used to maintain proper hip alignment by preventing the legs from turning inward, which can lead to hip dislocation or other complications.
- **Knee Extension:** The knee extender component is designed to support the knees in a slightly bent position, helping to reduce the pressure that lying flat can put on the lower back.

Model	Dim.	UOM
9913-24148H-P7	24"L x 14"D x 8.5"H / Pommel 7" x 3" x 3"	Each

Bedside Safety Mats

Meets CA 117 standards

Bedside Safety Mats, made from high-impact foam, help reduce the possibility of impact related injuries and provide full-length coverage of the bed exit area. Each safety mat is designed with a skid-resistant bottom, preventing slipping and aiding in ambulation, and features fire-resistant properties.

RECOMMENDED USE: Patients cognitively or physically disabled that are at risk of falls.

FALLSHIELD™ Bedside Safety Mat



The **NYOrtho FALLSHIELD™ Bedside Safety Mat** utilizes high-impact resistant, close-sealed foam that helps prevent patient injuries resulting from bedside falls. Aggressive non-slip bottom protects caregiver and patient from slipping and aids in ambulation. Bi-fold design allows for easy storage when not in use.

- Available in an easy to store Bi-Fold version.
- ¾ inch padding provides necessary impact attenuation performance requirements tolerance ratio.
- Impact absorbing durable molded foam/rubber prevents patient injury from falls.
- Low profile and beveled edge allow easy on and off for wheelchairs, beds and lifts.
- Ample size provides full coverage of the bed area.
- Non-slip textured top helps prevents caregiver/patient slipping and slide-resistant bottom keeps mat in place.
- Easy-to-clean fluid-proof cover.

Model	Style	Dim.	UOM
9576-3Q2470	Standard	¾x24x70 in.	Each
9576B-3Q2470	Bi-Fold	¾x24x70 in.	Each



Tri-Fold Beveled Mat W/Glow Strips



The **Tri-Fold Beveled Mat W/Glow Strips** provide full-length coverage of the bed exit area, reducing the possibility of impact related injuries. A tri-fold design makes it easy to store and a glow-in-the-dark beveled outer edge reduces the trip hazard and accommodates wheelchairs. Made from ULTRA high-density bonded foam, fire-resistant, and skid-resistant materials to provide safety, prevent slipping and aid in ambulation.

- One inch thick, ULTRA high-density bonded foam helps prevent patient injury from falls.
- Glow-in-the-dark strip illuminates the mat during the night to help prevent tripping.
- Low profile and beveled edge allow easy on and off for wheelchairs, beds and lifts
- Features a fluid-proof, antimicrobial, mildew-resistant, fire-resistant.
- Tri-fold design with handles easily slides under the bed when not in use.
- Skid-resistant bottom prevents slipping and aids in ambulation.
- Extra-durable 3-PLY vinyl, wipe-clean cover.



Model	Style	Dim.	UOM
9579TR-013070-GLO	Tri-Fold	1x30x70 in.	Each

Bedside Safety Mat



Made from high-impact foam, **NYOrtho Bedside Safety Mats** help reduce the possibility of impact related injuries, providing full-length coverage of the bed exit area. Bi-fold design easily slides under the bed when not in use. Designed with a skid-resistant bottom that prevents slipping and aids in ambulation.

- High-impact foam, available in 2, 3 and 4 inch, helps prevent patient injury from falls.
- Bi-fold design easily slides under the bed when not in use.
- Skid-resistant bottom prevents slipping.
- Extra-durable 3-PLY vinyl is antimicrobial, mildew-resistant and fluid-proof.
- Wipe-clean surface.



9579TR-1H2472

Model	Style	Dim.	UOM
9579-022460	Bi-Fold	2x24x60 in.	Each
9579-022466	Bi-Fold	2x24x66 in.	Each
9579-022472	Bi-Fold	2x24x72 in.	Each
9579-023666	Bi-Fold	2x36x66 in.	Each
9579-032466	Bi-Fold	3x24x66 in.	Each
9579-032472	Bi-Fold	3x24x72 in.	Each
9579-033072	Bi-Fold	3x30x72 in.	Each
9579-042460	Bi-Fold	4x24x60 in.	Each
9579-043668	Bi-Fold	4x36x68 in.	Each
9579TR-1H2472	Tri-Fold	1.5x24x72 in.	Each
Custom sizes available - any length, width or height.			

Bi-Fold Beveled Bedside Safety Mat



The **NYOrtho Bi-Fold Beveled Bedside Safety Mat** provides full-length coverage of the bed exit area, reducing the possibility of impact related injuries. Safety Mats promote easy storage with bi-fold design and feature a beveled outer edge to accommodate wheelchairs. Made from ULTRA high-density bonded foam, fire-resistant, and skid-resistant materials to provide safety, prevent slipping and aid in ambulation.

- One inch thick, ULTRA high-density bonded foam helps prevent patient injury from falls.
- Low profile and beveled edge allow easy on and off for wheelchairs, beds and lifts.
- Features a fluid-proof, antimicrobial, mildew-resistant, fire-resistant.
- Bi-fold design easily slides under the bed when not in use.
- Skid-resistant bottom prevents slipping and aids in ambulation.
- Extra-durable 3-PLY vinyl, wipe-clean cover.

Model	Style	Dim.	UOM
9579BV-012466	Bi-Fold	1x24x66 in.	Each

Low Bed Fall Safety Mat, Bi-Fold

The **NYOrtho Low Bed Fall Safety Mat** is designed for use in conjunction with a low bed to help reduce the possibility of impact related injuries, providing full-length coverage of the bed exit area. Bi-fold design for easy storage when not in use. Low beds are an integral part of fall prevention strategies in healthcare and home care settings, especially when used in conjunction with bedside mats that match the height of the bed.



- High-impact foam, available in 2, 3 and 4 inch, helps prevent patient injury from falls.
- Bi-fold design easily slides under the bed when not in use.
- Skid-resistant bottom prevents slipping.
- Extra-durable 3-PLY vinyl is antimicrobial, mildew-resistant and fluid-proof.
- Wipe-clean surface.

Model	Style	Dim.	UOM
9579-082470	Bi-Fold	8x24x70 in.	Each
9579-083076	Bi-Fold	8x30x76 in.	Each
9579-093076	Bi-Fold	9x30x76 in.	Each
9579-093660	Bi-Fold	9x36x60 in.	Each
9579-103076	Bi-Fold	10x30x76 in.	Each
9579-133076	Bi-Fold	13x30x76 in.	Each
9579-142468	Bi-Fold	14x24x68 in.	Each

Custom sizes available - any length, width or height.



Ultra Hip Protector

RECOMMENDED USE: Patients at risk of falling or hip fracture.

The **NYOrtho Ultra Hip Protector** combines our proven seamless fixation technology with energy-absorbing ultra-foam technology, to provide superior and comfortable protection for the trochanter.



- High impact-absorbing viscoelastic pads.
- 100% latex-free inner wear briefs are **ANTIMICROBIAL**, soft and comfortably worn under clothing.
- Low-profile pads are virtually undetectable.
- Standard sizes S-XXL.
- Unisex.

Model	Size	Fits Waist	Fits Hips	UOM
9962-S	Small	26-30 in.	33-37 in.	Each
9962-M	Medium	30-34 in.	37-41 in.	Each
9962-L	Large	34-38 in.	41-45 in.	Each
9962-XL	XL	38-42 in.	45-49 in.	Each
9962-2XL	2XL	42-46 in.	49-52 in.	Each

To insure proper fit, compare the waist measurement to hip measurement and select the size that is larger.

Hip Fractures

More than 95% of all hip fractures are caused by falls. The majority of hip fractures occur in adults over 65 years old. Fall-related injuries such as hip fractures are associated with increased mortality and a major cause of loss of functional capacity and impaired quality of life. Hip fractures are one of the largest health concerns in nursing homes as well as an extremely serious threat to the quality of life and the future health of elderly patients. Our Ultra Hip Protector offers distinctive advantages over competing products, through the use of unique materials.

THERAPY & AMBULATION

NYOrtho therapy products combine comfort, form and function to deliver lightweight and versatile designs.

Shoulder Orthosis Abductor

RECOMMENDED USE: Bankart repairs, SLAP lesion repairs, anterior dislocations, posterior capsule repairs, anterior repairs, rotator cuff repairs, and total shoulder reconstruction.



The **NYOrtho Shoulder Immobilizer** provides positional arm support and immobilization following upper extremity injuries and/or surgeries. This sling and pillow combination device can hold the shoulder/arm in a neutral position and help to prevent post-operative complications. New breathable spacer fabric and lightweight design enhances patient comfort.

- Easy to apply and to adjust arm position.
- Breathable padded spacer fabric for added comfort.
- Foam abduction pillow encased in breathable spacer fabric helps maintain neutral position.
- Comfortable low-profile foam pad cushions the sling against body when worn without the abductor pillow.
- Exercise ball for muscle and circulatory stimulation.
- Universal sizing – fits left and right arm.

Model	Size	UOM
8712-SM	Small/Medium	Each
8712-LXL	Large/XL	Each

Velpeau Arm Sling

RECOMMENDED USE: For non-surgical treatment of bone, muscle, or ligament injury to the hand, arm or shoulder requiring velpeau stabilization.



The **NYOrtho Velpeau Arm Sling** comfortably positions and stabilizes the arm to promote healing. Velpeau strap helps minimize shoulder rotation and provides comfortable support for either arm during recuperation from a broken bone, sprain, surgery, or for general arm support.

- Lightweight envelope sling features an open envelope design by the elbow and wrist.
- Velpeau arm sling helps control shoulder rotation with a stabilizing strap around the lower back.
- Easy to apply quick-release buckle.
- Ambidextrous, one size fits all – 7.5H x 16L in. envelope.
- Machine washable.
- Blue color may vary.

Model	Size	UOM
8700	Universal	Each

Arm Sling

RECOMMENDED USE: For non-surgical treatment of bone, muscle, or ligament injury to the hand, arm or shoulder.



The **NYOrtho Arm Sling** comfortably positions and supports the arm to promote healing. Provides comfortable support for either arm during recuperation from a broken bone, sprain, surgery, or for general arm support.

- Lightweight envelope sling supports arm.
- Available in an easy to apply hook and loop closure.
- Ambidextrous, sized S-XL.
- Machine washable.
- Blue color may vary.

Model	Size	Dim.	UOM
8701-S	Small	7.5x14 in.	Each
8701-M	Medium	7.5x16 in.	Each
8701-L	Large	7.5x18 in.	Each
8701-XL	XL	7.5x21 in.	Each

Gait Transfer Belts

Bedside Safety Mats, made from high-impact foam, help reduce the possibility of impact related injuries and provide full-length coverage of the bed exit area. Each safety mat is designed with a skid-resistant bottom, preventing slipping and aiding in ambulation, and features fire-resistant properties.

RECOMMENDED USE: Patients requiring assistance with transfer or ambulation.

Deluxe Gait Belt



The **NYOrtho Deluxe Gait Belt** offers the simplicity of a gait belt with the added support of easy to grasp hand loops. Provides caregivers with improved hold on the patient/resident to assist in safe transfer/ambulation and help prevent caregiver back injuries.

- Extra wide 4" high back for greater weight bearing leverage and comfort for the patient.
- Six handles provide attendants with a firm grip to assist in the safe transfer or ambulation of patients.
- Belt is secured with an adjustable quick-release buckle.
- Soft nylon is machine washable.
- Universal sizing.



Model	Fits	UOM
9507	28-55 in.	Each

Gait Belt



NYOrtho Gait / Transfer Belts provide caregivers with a leveraged hold on the patient/resident to assist in safe transfer/ambulation and help prevent caregiver back injuries.

- Nickel-plated, metal buckles offer superior hold and durability.
- Constructed of heavy-duty 2" cotton webbing and available in a variety of colors and patterns.
- Machine washable.

Style	Buckle	2x54 in.	2x60 in.	2x72 in.	UOM
Coral Reef	Metal	9503C-54	9503C-60	9503C-72	Each
Fiesta	Metal	9503F-54	9503F-60	9503F-72	Each
Lavender	Metal	9503L-54	9503L-60	9503L-72	Each
Natural	Metal	9503W-54	9503W-60	9503W-72	Each
Navy	Metal	9503N-54	9503N-60	9503N-72	Each
Pastel	Metal	9503P-54	9503P-60	9503P-72	Each
Stars & Stripes	Metal	N/A	9503S-60	N/A	Each
Natural	Delrin	9505W-54	9505W-60	9505W-72	Each
Pastel	Delrin	9505P-54	9505P-60	9505P-72	Each

Custom sizes available in **Natural** only - any length.



Palm Grips

RECOMMENDED USE: Individuals at risk of hand contractures.



NYOrtho Palm Grips provide stimulation and exercise for the hand ligaments and muscles, helping prevent hand contractures. Protect soft, irritated skin on fingers and palms. May be utilized as a diversionary device for individuals prone to pulling at IV or feeding tubes.

- Features a highly absorbent terry cloth fabric with **antimicrobial** protection that absorbs perspiration
- Latex-free elastic strap holds grip in place
- Available in two sizes: regular and large
- Machine washable



Model	Size	Dim.	UOM
9801-L	Large	6x3 in.	Each
9801-R	Regular	5x2.5 in.	Each

Fetal Monitoring Belts



NYOrtho Fetal Monitor Belts are used to keep fetal monitors in place. Made of soft, comfortable fabric, they are lightweight and cool for minimal perspiration.

- Soft, latex-free elastic belt comfortably holds fetal monitor in place.
- Available in button hole closure.
- Come in a set – 1 pink / 1 baby blue.

Model	Style	Dim.	UOM
9650	Button Hole	2.5x44 in.	Each

Soft Cervical Collar

RECOMMENDED USE: Mild neck pain, especially caused by a condition within the cervical spine.



NYOrtho Universal Cervical Collars helps provide relief from mild neck pain and stiffness. When used as prescribed by your physician, the NYOrtho Universal Cervical Collar should be worn for relief from pinched nerves or tightness in neck. It may also provide comfort after a whiplash injury.

- Breathable foam covered with soft, smooth knit material.
- Contoured shape for more comfortable fit.
- Available in 3" height.
- Adjustable velcro fastener promotes comfort.
- One size fits all.

Model	Height	UOM
8500	3 in.	Each

Double Lobe Cervical Pillow

RECOMMENDED USE: Helps provide comfort and support.



The **NYOrtho Double Lobe Cervical Pillow** is specially shaped to provide a uniform firmness while maintaining the proper curvature of the neck.

- Helps align the spine while promoting proper curvature of the neck.
- Medium density polyurethane foam.
- Cotton/poly cover.
- Machine washable.

Model	Dim.	UOM
9630	19x14 in.	Each

Post Hip Surgery Cushion

RECOMMENDED USE: Post hip surgery patients at risk of excessive flexion at hip joint.



The **NYOrtho Post Hip Surgery Cushion** is a lightweight and easy to use cushion that assists in the transition from sitting to standing from any chair. Fits most chairs, including straight backs and cushioned arm chairs.

- Super high density foam cushion provides necessary lift to reduce hip flexion when rising from seated position.
- Low shear water-resistant removable cover.
- Convenient carry strap.
- Meets CA 117
- Machine washable.

Model	Dim.	UOM
9590-151504	15x15x4 in.	Each
Custom sizes available in any length, width or height.		

Non-Surgical Reusable Isolation Gown Level 2 LB-P

RECOMMENDED USE: DO NOT use (I) in a surgical setting; (II) where significant exposure to liquid bodily or other hazardous fluids may be expected; or (III) in a clinical setting with high risk of infection.



Non-Surgical, Reusable Isolation Gowns are a cost conscious and environmentally responsible choice that provides protection to the wearer from the transfer of microorganisms and bodily fluids in low- or minimal-risk patient isolation situations and that are not intended for use during surgical procedures, invasive procedures, or when there is a medium or high risk of contamination.

- Meets the FDA and CDC guidelines for a Level 2 isolation gown.
- Polyester fabric tested to meet AAMI PB70 fluid barrier requirements for Level 2.
- Fabric maintains level 2 barrier up to 25 times commercial laundering.
- Overlap back with twill tie neck and back waist closures.
- Generously sized with longer sleeve length for comfortable arm movement with white ribbed cuffs.
- Reusable gown may be washed and dried – see washing instructions.
- Non-sterile.

Model	Size	UOM
FA-GWN-M2LB	Universal	Each
Custom sizes available in any length, width or height.		



PATIENT SAFETY

NYOrtho offers a variety of premium products to insure patient safety, helping protect against fall related injuries as well as restraints and restraint alternatives.

Door Guard

RECOMMENDED USE: For wandering patients; to deter unauthorized entry or exits.



The **NYOrtho Door Guard** is a visual reminder to deter wandering residents or patients from entering off limit zones and areas. The adjustable door guard provides fall management and cautionary discouragement for Alzheimer's and Dementia patients as well as for confused residents.

- Attaches easily with adhesive hook and loop strips to the door frame and removes easily in case of emergency.
- Alarm sounds when banner is detached.
- Battery operated with 9-Volt battery (included).
- Low battery indicator.
- Fits many size door frames.



Door Guard	Door Guard & Alarm	Dim.	UOM
9797-40	9797-40-ALRM	12x40 in.	Each
9797-48	9797-48-ALRM	12x48 in.	Each
9797-53	9797-53-ALRM	12x53 in.	Each
9797-60	9797-60-ALRM	12x60 in.	Each
	9797-ALRM		Each
Custom sizes available in any length or height.			

Padded Mitts

RECOMMENDED USE: Patients who are prone to pull at their IV line/catheter or who are prone to self-injury.



NYOrtho Padded Mitts protect hands from unwanted picking and scratching. Options include dual-sided padding or one-sided padding combined with breathable mesh fabric to provide additional comfort for hands. Wrist tie closure enables use without arm movement control straps.

- Available in lightweight mesh back that promotes air circulation and visual inspection (9370) or dual padded mitts (9371).
- Soft interior taslon fabric is gentle on the skin and water-resistant.
- Wrist ties allow for use without arm movement control straps.
- Half-inch padding cushions hand.
- No finger separators.
- One size fits all.
- Machine washable.



Model	Style	Size	UOM
9370	Padded / Mesh	Universal	Pair
9371	Dual Padded	Universal	Pair

Sheepskin Limb Holder

RECOMMENDED USE: Patients assessed to be at risk of interrupting treatment such as a skin condition, causes self-injury, or compromises wound site integrity.



The **NYOrtho Sheepskin Limb Holder** restraint comfortably secures wrists or ankles, limiting the motion of patients arms and legs in bed. Designed for at-risk patients, a non-slip buckle fastener provides a secure and easy to use closure. Designed with extra-long straps to allow fastening out of patient reach.

- Wrist and ankle restraint made of soft, high-pile, non-allergenic imitation sheepskin.
- Over the cuff adjustable slide allows wrist to rotate without chaffing.
- One size fits all.
- Machine washable.

Model	Size	UOM
9509-C	Universal	Pair

ORTHOPEDIC REHAB

NYOrtho offers a comprehensive line of orthopedic supports for a wide variety of applications. Supports help protect against muscle strain, improve lower back pain, and supplement abdominal strength.

Universal Rib Belt Hook Receptive

RECOMMENDED USE: Patients in need of stabilizing rib and sternum as a result of fractured or bruised ribs, intercostal muscle strains and costovertebral sprains.



The **NYOrtho Universal Rib Belt Hook Receptive** is specially designed to comfortably fit the female figure. Ideal for stabilizing rib and sternum fractures by limiting chest expansion through compression. Easy adjust-

- Hook compatible elastic permits a closure at any point on the circumference of the rib belt.
- Inner plush elastic is soft on the skin.
- Universal sizing.
- Latex-Free.



Model	Style	Fits Waist	UOM
8230-3045	Female	30-45 in.	Each
8230-4560	Female	45-60 in.	Each
8230-6075	Female	60-75 in.	Each
8232-3045	Male	30-45 in.	Each
8232-4560	Male	45-60 in.	Each
8232-6075	Male	60-75 in.	Each

Custom sizes available in any length or height.

Universal Abdominal Binder Hook Receptive

RECOMMENDED USE: Patients in need of postoperative support of the abdomen.



Whether recovering from surgery or injury, or just need extra support, **NYOrtho Universal Abdominal Binders Hook Receptive** can stabilize your core and improve muscle tone. Made of comfortable, flexible materials, abdominal binders secure post-surgical dressings and provide compression to the upper and lower abdomen while improving blood circulation and oxygen levels. This in turn reduces swelling and promotes healing.

- Hook compatible elastic permits a closure at any point on the circumference of the binder.
- Plush elastic is soft on the skin.
- Options: Available in 6", 9" or 12" height.
- Universal sizing. Bariatric Sizes available.
- Latex-Free.

Model 6 in.	Model 9 in.	Model 12 in.	Fits Waist	UOM
8102-3045	8100-3045	8101-3045	30-45 in.	Each
8102-4560	8100-4560	8101-4560	45-60 in.	Each
8102-6075	8100-6075	8101-6075	60-75 in.	Each
8102-7590	8100-7590	8101-7590	75-90 in.	Each

Custom sizes available in any length or height.



Tapered Abdominal Binder

RECOMMENDED USE: Patients in need of postoperative support of the abdomen.



Whether recovering from surgery or injury, or just need extra support, **NYOrtho Tapered Abdominal Binders** can stabilize your core and improve muscle tone. Made of comfortable, flexible materials, abdominal binders secure post-surgical dressings and provide compression to the upper and lower abdomen while improving blood circulation and oxygen levels. This in turn reduces swelling and promotes healing.

- One-panel plush elastic is soft on the skin.
- Tapered design for better fit.
- Effortless hook and loop closure.
- Available in 6", 9", 10", or 12" height.
- Standard sizing.
- Latex-Free.

Model 6 in.	Model 9 in.	Model 10 in.	Model 12 in.	Fits Waist	UOM
8000-S	8140-S	8142-S	8144-S	24-30 in.	Each
8000-M	8140-M	8142-M	8144-M	30-36 in.	Each
8000-L	8140-L	8142-L	8144-L	36-42 in.	Each
8000-XL	8140-XL	8142-XL	8144-XL	42-48 in.	Each
8000-2XL	8140-2XL	8142-2XL	8144-2XL	48-54 in.	Each
Custom sizes available in any length or height.					

Universal Abdominal Binder

RECOMMENDED USE: Patients in need of postoperative support of the abdomen.



Whether recovering from surgery or injury, or just need extra support, **NYOrtho Universal Abdominal Binders** can stabilize your core and improve muscle tone. Made of comfortable, flexible materials, abdominal binders secure post-surgical dressings and provide compression to the upper and lower abdomen while improving blood circulation and oxygen levels. This in turn reduces swelling and promotes healing.

- Soft medical foam panels for added comfort.
- Plush elastic is soft on the skin.
- Adjustable velcro fastener for comfortable fit.
- Options: Available in 3-panel 9", 1-panel 10" or 4-panel 12" height.
- Universal sizing.
- Latex-Free.

Model 9 in.	Model 10 in.	Model 12 in.	Fits Waist	UOM
8130-3045	8132-3045	8134-3045	30-45 in.	Each
8130-4662	8132-4662	8134-4662	46-62 in.	Each
8130-6075	8132-6075	8134-6075	60-75 in.	Each
8130-7590	8132-7590	8134-7590	75-90 in.	Each
Custom sizes available in any length or height.				

Lumbar Sacral Support DCDO

PROVIDES: Dual Compression Double Overlapping



1000/2000

Elastic body provides adjustable comprehensive support with semi-rigid optional removable pad or steel stays encased in a soft textile back.

- 10" high textile back tapering to 6" front elastic panels provides first layer of compression.
- Two double overlapping outer-belts provide additional compression to the abdominal area.
- Flexible polypro boning elements for added support.
- Double-ply textile back with soft flannel lining.
- Available with a removable foam pad or shaped steel stays.
- Adjustable velcro fastener for comfortable fit.
- Not made with natural rubber latex.
- Custom sizes available.

Foam Pad	Shaped Steel Stays	Size	Fits Waist	UOM
1000-S	2000-S	Small	26-30 in.	Each
1000-M	2000-M	Medium	30-34 in.	Each
1000-L	2000-L	Large	34-38 in.	Each
1000-XL	2000-XL	XL	38-42 in.	Each
1000-2XL	2000-2XL	2XL	42-46 in.	Each

Custom sizes available in any length or height.



1000



2000

Lumbar Sacral Support CC

PROVIDES: Criss-Cross



- 9" high all elastic back tapering to 4" front elastic panels.
- All elastic construction enfolds body contours comfortably.
- 4 flexible polypro boning elements for added support.
- Not made with natural rubber latex.

Model	Size	Fits Waist	UOM
3009-S	Small	26-30 in.	Each
3009-M	Medium	30-34 in.	Each
3009-L	Large	34-38 in.	Each
3009-XL	XL	38-42 in.	Each
3009-2XL	2XL	42-46 in.	Each

Custom sizes available in any length or height.



Universal Lumbar Sacral Support SCDS

PROVIDES: Single Compression Double Side-Pull



- 10" high elastic back tapering to velfoam panels
- Designed with double side pulls to provide additional compression to alleviate lower back pain.
- 4 flexible polypro boning elements for added support.
- Removable foam pad to relieve stress to the lower spine.
- Adjustable velcro fastener for comfortable fit.
- Not made with natural rubber latex.
- Universal sizing.

Model	Size	Fits Waist	UOM
5000-S	Small	21-26 in.	Each
5000-R	Regular	26-36 in.	Each
5000-XL	XL	36-46 in.	Each
5000-2XL	2XL	46-56 in.	Each

Custom sizes available in any length or height.



Lumbar Sacral Support DCS0

PROVIDES: Dual Compression Single Overlapping



Elastic body provides adjustable comprehensive support with semi-rigid optional removable steel stays encased in a elastic back.

- 10" high all elastic back tapering to 6" front elastic panels.
- Single overlapping belt provides additional abdominal support.
- All elastic construction enfolds body contours comfortably.
- 4 flexible polypro boning elements for added support.
- Available white or black.
- Not made with natural rubber latex.

White	Black	Size	Fits Waist	UOM
5500-S	5500BK-S	Small	26-30 in.	Each
5500-M	5500BK-M	Medium	30-34 in.	Each
5500-L	5500BK-L	Large	34-38 in.	Each
5500-XL	5500BK-XL	XL	38-42 in.	Each
5500-2XL	5500BK-2XL	2XL	42-46 in.	Each
Custom sizes available in any length or height.				



Universal Maternity Support

CAUTION: Please consult your physician before wearing this product.



The **NYOrtho Elastic Maternity Support** helps alleviate lower back pain due to pregnancy by supporting the abdominal and lumbar regions. Fully adjusts to accommodate abdominal growth throughout pregnancy. Elastic side panels add compression and ventilation. Hides well under clothing.

- Medical foam body applies over abdomen for optimum comfort.
- Plush elastic side panels for comfortable compression.
- Adjustable velcro fastener for comfortable fit.
- Universal sizing based on 6-14 dress size before pregnancy.
- Not made with natural rubber latex.

Model	Size	Fits Dress Size	UOM
8310-U	Universal	6-14 in.	Each
Custom sizes available in any length or height.			

OCCUPATIONAL & ERGONOMIC

NYOrtho Occupational Back Supports provide adjustable comprehensive support to help prevent the discomforts of low back pain. Helping stabilize lumbar, sacral, abdominal muscles and ligaments, relieving stress on the lower spine.

Breathable Spandex Back Belt



- 8" high breathable spandex back tapering to hook and loop front panels.
- 4" double outer-belt provides additional compression to abdominal area.
- 4 flexible polypro boning elements for added support.
- Interior gripper webbing prevents belt from sliding.
- Inner foam layer for extra comfort.
- 1½" removable adjustable suspenders.



Model	Size	Fits Waist	UOM
6027-S	Small	26-30 in.	Each
6027-M	Medium	30-34 in.	Each
6027-L	Large	34-38 in.	Each
6027-XL	XL	38-42 in.	Each
6027-2XL	2XL	42-46 in.	Each
Custom sizes available in any length or height.			

Deluxe Velfoam Back Belt



- 9" high all elastic back tapering to velfoam front-panels.
- 4" double outer-belt provides additional compression to abdominal area.
- 6 flexible polypro boning elements for added support.
- Inner foam layer for extra comfort.
- 1½" adjustable non-stretch suspenders.



Model	Size	Fits Waist	UOM
7100-S	Small	26-30 in.	Each
7100-M	Medium	30-34 in.	Each
7100-L	Large	34-38 in.	Each
7100-XL	XL	38-42 in.	Each
7100-2XL	2XL	42-46 in.	Each
Custom sizes available in any length or height.			

Deluxe Ventilated Elastic Back Belt



- Breathable ventilated elastic for comfort.
- 9" high elastic back tapering to 6" elastic front-panels.
- Overlapping 4" single outer-belt provides additional compression to abdominal area (Low-Profile 3").
- 4 flexible polypro boning elements for added support.
- Interior gripper webbing prevents belt from sliding.
- 1½" adjustable suspenders.



Model	Size	Fits Waist	UOM
7505-S	Small	26-30 in.	Each
7505-M	Medium	30-34 in.	Each
7505-L	Large	34-38 in.	Each
7505-XL	XL	38-42 in.	Each
7505-2XL	2XL	42-46 in.	Each
Custom sizes available in any length or height.			

Deluxe Ventilated RF Elastic Back Belt



- Breathable ventilated elastic for comfort.
- 9" high elastic back tapering to 6" elastic front-panels.
- Overlapping orange reflective 4" single outer-belt provides additional compression to abdominal area (Low-Profile 3").
- 4 flexible polypro boning elements for added support.
- Interior gripper webbing prevents belt from sliding.
- 1½" orange reflective adjustable suspenders.

Model	Size	Fits Waist	UOM
7505RFL3-S	Small	26-30 in.	Each
7505RFL3-M	Medium	30-34 in.	Each
7505RFL3-L	Large	34-38 in.	Each
7505RFL3-XL	XL	38-42 in.	Each
7505RFL3-2XL	2XL	42-46 in.	Each
Custom sizes available in any length or height.			

Deluxe Breathable Spandex Back Belt



- 9" high breathable spandex back tapering to hook and loop front panels.
- Vinyl reinforced side-pulls.
- 4" double outer-belt provides additional compression to abdominal area.
- 4 flexible polypro boning elements for added support.
- Interior gripper webbing prevents belt from sliding.
- 1½" adjustable suspenders.

Model	Size	Fits Waist	UOM
6800BK-S	Small	26-30 in.	Each
6800BK-M	Medium	30-34 in.	Each
6800BK-L	Large	34-38 in.	Each
6800BK-XL	XL	38-42 in.	Each
6800BK-2XL	2XL	42-46 in.	Each
Custom sizes available in any length or height.			

Elastic Wrist Wrap



9450-BK

The NYOrtho Wrist Wrap protects weak or injured wrists by providing reliable support and limiting wrist motion. Designed for all-day comfort. Fits discreetly under clothing.

- Durable 3" elastic for support and comfort
- Thumb loop for ease of application
- Adjustable strap for custom fit and level of support
- Available without thumb loop (9460)
- Size: Universal / Ambidextrous



9460-23

Model	Style	Size	UOM
9450-BK	W/Thumb Loop	Universal	Each
9460-23	Wrap	Universal	Each

PROTECTIVE SAFETY WEAR

Workplace injuries are costly to both employer and employee. NYOrtho's protective and safety wear products ensure workplace injuries don't result, increasing health and efficiency in the work force.



ANSI 107-2015 CLASS 2 SAFETY VEST **American National Standard for High-Visibility Safety Apparel** **TYPE R SAFETY VEST** Roadway and Temporary Traffic Control Zones



HI-VIS RETRO REFLECTIVE MATERIAL

At least 201 in²

Minimum Width: 1.38" / Actual: 2"

360° COVERAGE

One Horizontal Band &
Two Vertical Bands

DURABLE, BREATHABLE POLYESTER MESH FABRIC

Hi-Vis Background Material
At least 775 in²

HOOK & LOOP CLOSURE

Easy-on, Easy-off

ADJUSTABLE SIDE STRAPS

Expandable side straps in two sizes

SF-V2000: M/L 32"-44" size chest

SF-V2001: L/XL 44"-58" size chest



CRYOGEAR Hi-Visibility Cold Weather Jacket

The **NYOrtho Cryogear Cold Weather Jacket** is a thoughtfully designed insulated winter parka jacket that combines durable polyester fabric with a modern fit and stylish look. Ideal for cold weather conditions, this jacket provides excellent thermal insulation with heavyweight lining. ANSI 107-2015 Class 3 hi-visibility standard keeps you safe as well as warm.

- Hi-visibility, water-resistant outer shell keeps you dry and safe.
- Thermal insulation provides superlative warmth and comfort.
- Heat reflective lining reflects heat back into the body ensuring maximum warmth and comfort.
- ID pocket for easy identification.
- Hook & loop adjustable width cuff.
- 4 zippered pockets provide ample storage.
- Internal zippered pocket provides safe and secure storage.

Model	Size	UOM
SF-CT-1001	Small	Each
SF-CT-1002	Medium	Each
SF-CT-1003	Large	Each
SF-CT-1004	XL	Each
SF-CT-1005	2XL	Each
SF-CT-1006	3XL	Each
SF-CT-1007	4XL	Each
SF-CT-1008	5XL	Each

ANSI 107-2015 CLASS 3 HI-VISIBILITY JACKET

American National Standard for High-Visibility Safety Apparel

TYPE R HI-VISIBILITY COAT

Roadway and Temporary Traffic Control Zones

CRYOGEAR Hi-Visibility Phase Vest

The **NYOrtho Cryogear Hot Weather Vest** utilizes phase technology to keep the wearer cool and comfortable when working conditions get hot. Four specially formulated cool-packs maintain a comfortable 58° F for up to 3 hours and only require 45 minutes to freeze.. Rotate a second set of cool packs for continuous relief. Universal sizing eliminates need for multiple sizes with an adjustable side strap. ANSI 107-2015 Class 2 hi-visibility standard keeps you safe as well as cool.

Model	Size	Fits Chest
SF-PH-2000	M/L	32-44 in.
SF-PH-2001	L/XL	44-58 in.

Picture unavailable.

Protective Knee Pads



NYOrtho Protective Knee Pads offer plenty of protection on hard, rough and even rocky surfaces. Comfortable foam padding and a tough fabric cover make these knee pads ideal for all day projects such as laying carpet or setting tile. The tough polypropylene caps are riveted in place and the knee pads secure in place with easy hook and loop straps.

- Hard polypro caps are riveted in place for protection on hard, rough, even rocky surfaces.
- Comfortable foam padding.
- Adjustable hook and loop closures.
- Durable polyester fabric cover.

Model	Size	UOM
SF-KP1	Universal	Each

OTG Safety Glasses



- Extremely comfortable, non slip design.
- Fits securely over prescription glasses.
- Anti-fog, anti-scratch lens.
- 99.9% UV protection.

Model	Size	UOM
SF-OTG-1000	Universal	Each

Cut Resistant Gloves

RECOMMENDED USE: Light/Medium Cut Hazards - material handling, small parts assembly with sharp edges, packaging, warehouse, general purpose, forestry, construction, pulp & paper, automotive assembly.



NYOrtho Cut Resistant Gloves are manufactured from a blended composite fiberglass, providing one of the highest performance cut resistant gloves on the market. Covered in a durable and lightweight HPPE fiber and wrapped in nylon/lycra, NYOrtho Cut Resistant Gloves are comfortable to wear, yet durable. Each size has a unique color band for easy size identification.

- Cut resistant and ergonomic design protects hand while maintaining comfort.
- Color-coded elasticized edge reduces fraying and provides a visual aid for sizing.
- Gloves can be laundered to increase wear time.
- EN388 Level 5 Cut Resistance / ANSI Cut Level A2.

Model	Size	UOM
SF-CRG-S	Small	Each
SF-CRG-M	Medium	Each
SF-CRG-L	Large	Each
SF-CRG-XL	XL	Each



63 FLUSHING AVENUE
BUILDING 280, SUITE 702
BROOKLYN, NEW YORK 11205

800-699-2358
www.nyorthousa.com

NEW PRODUCT DIRECTORY



ZERO-G Heel Cloud
pages 6-7



Amputee Gel-Foam Cushion
page 21



WaveForm Gel-Foam Cushion
page 26



TRIFORM Ultra
page 27



Geri-Chair Gel-Foam Comfort Seat
page 28



BalancePro Lateral Support
page 32



CurveFit Wheelchair Support
page 33



Foot & Calf Pad
page 33



CushionGuard Wheelchair Footrest Sleeve
page 37



Lateral Position Knee Abductor
page 41



Hip Abductor & Knee Extender Combo
page 43



Tri-Fold Beveled Mat W/Glow Strips
page 44



Low Bed Fall Safety Mat
page 46

NYOrtho

CUSTOMIZING COMFORT
THROUGH DESIGN & INNOVATION